→REF KST11003TP  
KST11003GP  
KST11003BW

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→REF KST11003TP  
KST11003GP  
KST11003BW

### Intended Use

The PreventID® CalDetect® 50/200 is an immunological rapid test for the determination of calprotectin in faeces.

### Introduction

Calprotectin (MRP 8/14) is a heteromer of two calcium-binding proteins (MRP8 and MRP 14) present in the cytoplasm of neutrophils and expressed by the membranes of monocytes. It plays a central role in neutrophil defense and may be detected in serum, body fluids or stool as inflammatory marker. Faecal calprotectin has been established as marker of inflammatory bowel diseases (IBD). It allows a reliable differentiation between IBD and functional intestinal diseases (e. g. irritable bowel syndrome, IBS). In addition, calprotectin is ideal for monitoring disease activity (e. g. of M. Crohn) and for early relapse detection.

The two cut-offs at 50 and 200 µg/g of the PreventID® CalDetect® 50/200 test enable not only the differentiation between inflammatory (≥50 µg/g) and non-inflammatory (<50 µg/g) diseases but also the individual therapy control of patients with inflammatory bowel diseases (IBD). A faecal calprotectin concentration above 200 µg/g in drug-treated IBD patients indicates a relapse and can be utilized to adjust therapy accordingly. A periodic application of the PreventID® CalDetect® 50/200 in IBD patients is therefore an effective tool for individual therapy control.

### Test Principle

PreventID® CalDetect® 50/200 is an immunological lateral flow test, which detects human calprotectin via gold-conjugated anti-calprotectin antibodies. A control line shows the correct test procedure. Two test lines allow calprotectin to be detected in the sample.

### Materials

#### Materials Provided

1. Test devices, individually packed **TEST**
2. Sample collection tubes with extraction buffer solution and sample collection stick **TUBE**
3. Paper stool catcher for faecal samples
4. Instruction sheets for sample collection for the patients
5. Manual

**Materials Required but not Provided:** Timer or stop watch

### Storage and Stability

The PreventID® CalDetect® 50/200 test should be stored at room temperature. The rapid test device is susceptible to humidity and high temperatures. Therefore, this test should be protected from extreme temperature and should be run immediately after the opening of the pouch. This test should not be run after the expiry date.

### Precautions

1. For *in vitro* diagnostic use only.
2. Do not eat or smoke while handling specimen.
3. Wear protective gloves and wash hands thoroughly after performing the test.
4. Avoid splashing or aerosol formation while handling specimen and performing the test.
5. All samples and materials used should be treated as potentially infectious and disposed in a biohazard container. Clean all contaminated objects and surfaces carefully.
6. Do not use test if the pouch is torn or if the membrane of the test device is visibly damaged. Note the expiry date.
7. Read the instruction carefully before performing the test.
8. Do not mix reagents from different lots.
9. If you have any questions please contact Preventis GmbH.

### Specimen Collection

1. Collect a stool sample with the paper stool catcher attached to your toilet according to the manual. The sample must not get in contact with water or urine.
2. Unscrew the cap of the sample collection tube and stick the attached sample collection stick **in one go at three different sites into the faeces**. Only the amount of stool that sticks to the grooves of the sample collection stick should be transferred to the sample collection tube.
3. Now retract the sample collection stick with the adhering faecal sample and insert it **only once** into the sample collection tube containing an extraction buffer solution.  
**Please note:** A repeated transfer of stool into the sample collection tube compromises the test performance!
4. Screw cap on firmly and shake well. This defined stool sample solution is now ready to use for the test.
5. If the test is not run on the day of sample collection, the sample collection tube should be stored at 2–8 °C, but no longer than 7 days.

### Test Procedure

1. Remove the test device from the pouch and place it on a flat dry surface. The oval sample application window at the one end of the test device should be at the right side (Fig. 1). The test device should immediately be used after opening the pouch.


Fig. 1: PreventID® CalDetect® 50/200 test device

2. If necessary bring sample collection tube to room temperature after sample collection and shake again.
3. Break off the tip of the sample collection tube carefully (avoid dripping). Squeeze **3 drops** of the extracted sample into the oval sample application window on the right side of the test device (by gently pressing the sample collection tube in the middle).
4. In a properly working test, a violet line will pass through the square result window in the middle of the test device.
5. The result should be interpreted **10 minutes** after the last drop has been placed.

### Test Interpretation (Fig. 2)

A solitary red control line (C) in the result window indicates that the test has run correctly. Depending on the concentration of calprotectin in the sample, test lines (T1, T2) will appear to the right of C.

**Negative:** Only the red control line (C) is visible. The test has run correctly, an intestinal inflammation could not be detected.

**Positive: Calprotectin concentration ≥ 50 µg/g:** The control line (C) and the test line T1 are visible. An inflammatory process has been detected.

**Calprotectin concentration ≥ 200 µg/g:** The control line (C) and both test lines T1 and T2 are visible. A high-grade inflammatory process has been detected.

**Invalid:** The test result is invalid if the control line (C) does not appear (even if test lines are visible).


Fig. 2: PreventID® CalDetect® 50/200 test interpretation

### Test Limitations

Although the PreventID® CalDetect® 50/200 is very accurate in detecting calprotectin a low incidence of false results may occur. Other clinically available tests are required if questionable results are obtained.

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by a physician after all clinical and laboratory findings have been evaluated.

### References:

1. Puolanne AM et al. (2016), Eur J Clin Invest 46(10):825-32
2. Aschauer GJM et al. (2010), Labmed 2010, P054
3. Tursi A et al. (2009), Int J Colorectal Dis 24:49–55

### Short instructions

#### PreventID® CalDetect® 50/200

1. Collect a faecal sample with the aid of the sample collection tube and the sample collection stick as described in the instruction.
2. Shake the solution in the sample collection tube very thoroughly. Unpack the test device.
3. Break off the tip of the sample collection tube carefully.
4. Squeeze **3 drops** of the extracted sample into the oval sample application window.
5. Interpret the test **after 10 minutes**.

Status: 2019-11-14

US: all products: Research Use Only. Not for use in diagnostic procedures.

PreventID® and CalDetect® are trademarks of Immundiagnostik AG, Bensheim.

Temperature limitation	Manufacturer
In vitro diagnostic device	Lot number
Catalogue number	Expiry date
To be used with	Do not reuse
Read user instructions	Contains sufficient for <n> tests
Keep away from sunlight	

Distributed by:  
**Preventis GmbH**  
Stubenwald-Allee 8a  
64625 Bensheim, Germany  
Phone: +49 6251 70711-0  
Fax: +49 6251 70711-299  
info@preventis.com  
www.preventis.com