→REF KST11004TP
KST11004GP
KST11004BW

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→REF KST11004TP
KST11004GP
KST11004BW

Intended Use

The PreventID® CalScreen® is an immunological rapid test for the detection of calprotectin in faeces. The determination of faecal calprotectin allows the differentiation between organic intestinal diseases (e.g. inflammatory bowel syndrome) and functional intestinal diseases (e.g. inflammatory bowel diseases).

Introduction

Calprotectin (MRP 8/14) is present in the cytoplasm of neutrophils and plays a central role in immune defense. Upon neutrophil activation or endothelial adhesion of monocytes, calprotectin is released and may be detected in serum, body fluids or stool as a potentially useful clinical inflammatory marker. The acute phase protein shows a high stability in faeces (stable for one week at room temperature) and has been established as a faecal marker of inflammatory bowel diseases (IBD).

Calprotectin enables a reliable differentiation between organic intestinal diseases (e.g. chronic inflammatory diseases, infectious diseases) and functional intestinal diseases (e.g. irritable bowel syndrome). Furthermore, the parameter is used to monitor disease activity in IBD (e.g. Crohn's disease or ulcerative colitis). Calprotectin serves as a positive predictive marker for invasive pathogens in screenings for infectious diarrhea.

The use of this test with a cut-off at 50 µg/g calprotectin enables the rapid screening of patients with chronic intestinal problems and the differentiation in inflammatory (≥50 µg/g) and non-inflammatory (<50 µg/g) diseases.

Test Principle

The PreventID® CalScreen® is an immunological lateral-flow test which detects human calprotectin via gold-conjugated anti-calprotectin antibodies. A control line indicates a correct test run, a test line provides evidence for the presence of calprotectin in the sample at a concentration of 50 µg/g or more.

Materials

Materials Provided

1. Test devices, individually packed **TEST**
2. Sample collection tubes with extraction buffer solution and sample collection stick **TUBE**
3. Paper stool catcher for faecal samples
4. Instruction sheets for sample collection for the patients
5. Manual

Materials Required but not Provided: Timer or stop watch

Storage and Stability

The PreventID® CalScreen® test should be stored at room temperature. The rapid test device is susceptible to humidity and high temperatures. Therefore, this test should be protected from extreme temperature and should be run immediately after the opening of the pouch. This test should not be run after the expiry date.

Precautions

1. For *in vitro* diagnostic use only.
2. Do not eat or smoke while handling specimen.
3. Wear protective gloves and wash hands thoroughly after performing the test.
4. Avoid splashing or aerosol formation while handling specimen and performing the test.
5. All samples and materials used should be treated as potentially infectious and disposed in a biohazard container. Clean all contaminated objects and surfaces carefully.
6. Do not use test if the pouch is torn or if the membrane of the test device is visibly damaged. Note the expiry date.
7. Read the instruction carefully before performing the test.
8. Do not mix reagents from different lots.
9. If you have any questions please contact Preventis GmbH.

Specimen Collection

1. Collect a stool sample with the paper stool catcher attached to your toilet according to the manual. The sample must not get in contact with water or urine.
2. Unscrew the cap of the sample collection device and stick the attached sample collection stick **in one go at three different sites into the faeces**. Only the amount of stool that sticks to the grooves of the sample collection stick should be transferred to the sample collection device.
3. Now retract the sample collection stick with the adhering faecal sample and insert it **only once** into the sample collection device containing an extraction buffer solution.
Please note: A repeated transfer of stool into the sample collection device compromises the test performance!
4. Screw cap on firmly and shake well. This defined stool sample solution is now ready to use for the test.
5. If the test is not run within one day of sample collection, the sample collection device should be stored at 2–8°C, but not longer than 7 days.

Test Procedure

1. Remove the test device from the pouch and place it on a flat dry surface. The oval sample application window at the one end of the test device should be at the right side (Fig. 1). Use test device immediately.


Fig. 1: PreventID® CalScreen® Test device

2. If necessary, bring sample collection device to room temperature after sample collection and shake again.
3. Break off the tip of the sample collection tube carefully (avoid dripping). Squeeze **3 drops** of the extracted sample into the oval sample application window on the right side of the test device (by gently pressing the sample collection tube in the middle).
4. In a properly working test, a violet line will pass through the square result window in the middle of the test device.
5. The result should be interpreted **10 minutes** after the last drop has been placed.

Test interpretation (Fig. 2)

A solitary red control line (C) in the result window indicates that the test has run correctly. Depending on the calprotectin concentration, a test line (T) will appear to the right of the control band (Fig. 2a-c).

Negative:


Only the red control line (C) appears (Fig. 2a). The test has run correctly, no intestinal inflammation has been detected.

Positive:

Calprotectin concentration ≥50 µg/g: Control line (C) and test line (T) are visible (Fig. 2b). An intestinal inflammation has been detected.

Invalid:

The test is invalid if no control line (C) appears, even if a test line (T) is visible (Fig. 2c).

Test Limitations

Although the PreventID® CalScreen® is very accurate in detecting calprotectin a low incidence of false results can occur. Other clinically available tests are required if questionable results are obtained.

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

References:

1. Puolanne AM et al. (2016), Eur J Clin Invest 46(10):825-32
2. Aschauer GJM et al. (2010), Labmed 2010, P054
3. Tursi A et al. (2009), Int J Colorectal Dis 24:49–55

Short instructions PreventID® CalScreen®

1. Collect a faecal sample with the aid of the sample collection tube and the sample collection stick as described in the instruction.
2. Shake the solution in the sample collection tube very thoroughly. Unpack the test device.
3. Break off the tip of the sample collection tube carefully.
4. Squeeze **3 drops** of the extracted sample into the oval sample application window.
5. Interpret the test **after 10 minutes**.

Status: 2018-11-29

US: all products: Research Use Only. Not for use in diagnostic procedures.

PreventID® and CalScreen® are trademarks of Immundiagnostik AG, Bensheim.

Temperature limitation	Manufacturer
In vitro diagnostic device	Lot number
Catalogue number	Expiry date
To be used with	Do not reuse
Read user instructions	Contains sufficient for <n> tests
Keep away from sunlight	

Distributed by:

Preventis GmbH
Stubenwald-Allee 8a
64625 Bensheim, Germany
Phone: +49 6251 70711-0
Fax: +49 6251 70711-25
info@preventis.com
www.preventis.com