# PreventID® CC

## Testanleitung

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## PreventID® CC

### Test instructions

PreventID® CC rapid test for the determination of occult blood in stool – for self-testing

### Prevention and early detection of colon cancer

#### Intended use

The PreventID® CC is an easy-to-use immunological rapid test for the qualitative detection of human faecal blood in stool, intended for self-testing.

#### Introduction

In general, colon cancer develops over a period of years from benign polyps in the intestinal wall. Not every polyp progresses to a malignant tumour but nevertheless should be removed as a precaution.

Polyps tend to bleed temporarily and hence can be discovered at an early stage without invasive surgery merely by detecting blood traces in a stool sample.

Conventional guaiac tests lack the sensitivity and specificity to reliably detect occult blood. Furthermore, these tests require a strict diet before the procedure.

The PreventID® CC test is able to detect occult blood in far lower concentrations compared to the conventional guaiac test. Moreover, the specificity of the PreventID® CC test is not influenced by interfering food components.

Finally, the PreventID® CC test result is easier to interpret than the guaiac tests result.

#### Test principle

This test is based on a sandwich immunoassay whose unique combination of monoclonal antibodies

permits the selective detection of occult blood in faecal samples. By this method, human occult blood can be detected with high sensitivity and specificity.

Gold-labelled anti-hemoglobin antibodies bind to hemoglobin (if present) in the sample. These complexes migrate over the test membrane until they pass the test zone which consists of immobilized anti-hemoglobin antibodies. These secondary antibodies bind the gold-labelled antibody-hemoglobin complexes and thus generate a coloured (test) line (positive result).

In order to prevent false-negative results, an internal (control) line consisting of anti-IgG antibodies has been added. These bind excess gold-labelled antibodies and retain them on the control line. A coloured control line thus indicates that the test has run properly and the result is valid.

#### Test characteristics

**Analytical Sensitivity:** 50 ng/mL. The test indicates the result as positive with occult blood higher than 50 ng/mL.

**Analytical Specificity:** sample were spiked with the following organisms and substances (e. g. bacteria, yeasts, drugs, food, food additives), tested with PreventID® CC and showed no cross-reactivities or interferences: Acinetobacter calcoaceticus, Proteus vulgaris, Salmonella typhi, Acinetobacter

spp, Staphylococcus aureus, Candida albicans, Neisseria gonorrhoeae, Escherichia coli, Neisseria catarrhalis, Gardnerella vaginalis, Neisseria meningitidis, Streptococcus faecalis, Neisseria lactamica, Streptococcus faecium, Pseudomonas aeruginosa, Trichomonas vaginalis; Pseudoglobins from beef, fish, chicken, horse or sheep; broccoli, cantaloupe, horseradish, cauliflower, vitamin C, iron; Acetaminophen, Acetyl salicylic acid, Ascorbic acid, Atropine, bilirubin, caffeine, creatinine, gentesic acid, glucose, ketones, Mestranol, nitrite, Penicillin, sodium- and lithium heparin.

Only with hemoglobin concentrations higher than 100,000 ng/mL the test shows false negative results. PreventID® CC has a **clinical sensitivity** of 75 % (positive samples are evaluated positive) and a **clinical specificity** of 99 % (negative samples are evaluated as negative) relative to colonoscopy.

#### Materials provided

- 1 test device [TEST] (with drying agent; not required for the test)
- 1 sample collection tube with buffer and sample collection stick [TUBE]
- 1 paper stool catcher for faecal samples
- test instructions

**Additionally required:** Timer or stop watch.

#### Storage, Stability and Disposal

The PreventID® CC test should be stored at room temperature (4–30 °C). Do not freeze. The test device is susceptible to humidity and high temperatures. Therefore, this test should be protected from extreme temperature and should be run immediately after opening of the pouch. This test should not be run after the expiry date. Dispose all test components and samples in conventional garbage, the paper stool catcher into the toilet.

#### Precautions

1. For *in vitro* diagnostic use only. Only for application outside the body.
2. Read the test instructions before using the test.
3. Do not eat or smoke while handling specimen. Wear protective gloves if available and wash hands thoroughly after performing the test.
4. Avoid splashing or aerosol formation.
5. Clean all contaminated surfaces carefully.
6. Do not use test if the pouch is torn or if the membrane of the rapid test device is visibly damaged. Do not use the test after expiry date.
7. If your skin or eyes had contact with the buffer immediately wash with plenty of water.
7. Keep away from children.
8. If you have questions please contact Preventis GmbH.

#### Sample preparation

1. Sample collection should not be performed during or within three days of a menstrual period, or if you suffer from bleeding gums or haemorrhoids or blood in the urine. These conditions may cause false-positive results.
2. No dietary requirements.
3. Alcohol, aspirin and some medications sometimes cause gastrointestinal irritation or bleeding. Consult your doctor prior to the test whether and to which extent you should temporarily reduce the intake of such substances.

#### Sample collection

1. Collect a faecal sample with the paper stool catcher attached to your toilet according to the manual. Please note: The test result can be adversely affected when faeces come in contact with toilet water, urine or toilet bowl freshener.
2. Unscrew the cap of the sample collection tube and stick the attached sample collection stick **in one go at three different sites into the faeces** (Fig. 1). Only the amount of stool that sticks to the grooves of the sample collection stick should be transferred to the sample collection tube.
3. Now retract the sample collection stick with the adhering faecal sample and insert it **only once** into the sample collection tube containing an

extraction buffer solution.

**Please note:** A repeated transfer of stool into the sample collection tube compromises the test performance!

4. Screw cap on firmly and shake well. This defined stool sample solution is now ready to use for the test.
5. The paper stool catcher can be flushed down the toilet.

#### Test procedure

After specimen collection bring all test components to room temperature.

1. Remove test from its pouch and place on a flat dry surface. The oval sample opening at the one end of the test device should be at the right side (Fig. 2).
2. Break off the tip of the sample collection tube carefully as indicated in Fig. 2 A (avoid dripping by using a soft paper/toilet paper and holding the tip upwards). Squeeze the tube and apply **two drops** of the extracted sample into the sample opening on the right end of the test.
3. In a properly working test, a violet line will pass through the square result window in the middle of the test device.
4. The result should be interpreted **at 10 minutes**.

**Beware:** The above interpretation time is based on reading the test results at 15–30 °C.

#### Test interpretation

1. A coloured control (C) line at the left side of the result window indicates that the test has run correctly.
2. The test result is indicated by the appearance or non-appearance of a coloured test (T) line on the right side of the result window.

**Positive Result:** If both test (T) and control (C) lines appear in the result window, the test is positive, even if the test line is weak and regardless of which line is generated first (Fig 2B). Human occult blood has been detected in the stool sample. This does not automatically indicate the presence of colon cancer! The bleeding cause however should be investigated.

**Negative Result:** A solitary control (C) line in the result window indicates a negative result (Fig. 2 C). No human occult blood had been detected in the stool sample.

**Invalid Result:** If no control (C) line is visible after performing the test, the test result is considered invalid. The test result is also considered invalid if only the test (T), but no control (C) line appears (Fig. 2 D). In this case, the test must be repeated with a new test device.

**We recommend to consult a doctor promptly for a differential diagnosis following a positive result.** Please make no medically-relevant decision without consulting your attending physician.


Fig. 2: Test interpretation

#### Limitations of the test

- The origin of blood of faeces may be other than from colorectal bleeding such as haemorrhoids, blood in urine or gastric irritations.
- Negative results do not rule out bleeding because some polyps or adenomas bleed intermittently and blood is not distributed homogeneously within the faecal sample and may thus be missed during the sample collection.
- Colorectal polyps may not bleed during very early stages and sometimes bleed intermittently. It is recommended to check annually for occult blood in stool.
- Other clinical tests are indicated in the case of unspecific results. As with other diagnostic tests, a definitive clinical diagnosis should not be based on a single test result, but should be made by a qualified professional after evaluation of all clinical and laboratory findings.

#### Literature:

- Hoepfner N et al. (2006) Comparative evaluation of a new bedside faecal occult blood test in a prospective multicentre study. Aliment Pharmacol Ther 23: 145-154
- Schröder O et al. (2005) Klinische Validierung eines neuen immunologischen Bedside-Tests auf okkultes Blut im Stuhl. Z Gastroenterol 43:947 (P509)

2018-07-17

 Temperature limitation	 Manufacturer
 In vitro diagnostic device	 Lot number
 Catalogue number	 Expiry date
 Keep away from sunlight	 Do not reuse
 Read user instructions	 Contains sufficient for <n> tests

Learn more about colon cancer prevention with PreventID® CC at [www.prevent-colon-cancer.com](http://www.prevent-colon-cancer.com)