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EN Manual

Intended use

PreventID® Dermatophyte is an immunochromatographic rapid test for the qualitative determination of dermatophyte-derived antigens in nails.

Introduction

A fungal infection of nails, also called nail mycosis or onychomycosis, is an infection of the toenails or fingernails, usually caused by dermatophytes (filamentous fungi). Internal (endogenous) disposition factors could be for example metabolic diseases, genetic predispositions causing a susceptibility to onychomycosis, and immune deficiencies. Especially patients with diabetes mellitus, circulatory disorders, nail eczema and nail psoriasis are susceptible for fungal infections. Favouring is a moist-warm environment, toenails, therefore, are more frequently affected than fingernails. Another cause might also be a prolonged therapy with antibiotics.

Test Principle

PreventID® Dermatophyte is based on an immunochromatographic method for the determination of dermatophyte-antigens in nails infected with dermatophytes. The test device contains a fixed antibody against dermatophyte antigen in the test zone and a gold-conjugated monoclonal antibody against dermatophyte-antigen in the reagent zone. With the extraction solution the dermatophyte antigen is solubilised from the nail substance before the test run. In case there is dermatophyte-antigen in the extraction solution this antigen binds during the test run to the colloidal-gold antibody. This complex migrates to the fixed second monoclonal antibody forming a coloured precipitate at the test line and indicating a positive result – an infection with dermatophytes. A coloured control line is formed as internal control and proves the correct test run.

Material

- 5 / 10 test strips **[TEST]**, single packed*
- 1 extraction buffer **[BUF]**
- 5 / 10 test tubes **[TUBE]**
- 5 / 10 stir rods
- 1 tube rack
- 1 manual

*To prevent the test strip from kinking, the test pack contains a stabilising cardboard strip that can be disposed of together with the packaging after having taken out the test strip.

Material required but not provided: Timer or stop watch, gloves, clippers or scissors

Storage and Stability

Store the test between 2°C and 30°C; do not freeze. The test device is sensitive to humidity, direct sunlight as well as to heat. Perform the test immediately after removing the test device from the pouch. Do not use it beyond the expiry date.

Precautions

If the extraction buffer comes in contact with the eyes, mouth, or skin, rinse thoroughly with running water as first aid, and seek medical treatment if necessary.


- For *in vitro* diagnostic use only.
- Do not eat or smoke while handling specimen.
- Wear protective gloves and wash hands thoroughly after performing the test.
- Avoid splashing or aerosol formation while handling specimen and performing the test.
- All samples and materials used should be treated as potentially infectious and disposed in a biohazard container. Clean all contaminated objects and surfaces carefully.
- Do not use test if the pouch is torn or if the membrane of the test device is visibly damaged.
- Read the instruction carefully before performing the test.
- Do not mix reagents from different lots.
- If you have any questions please contact **Preventis GmbH**.

Specimen Preparation and Specimen Collection

This kit is intended for the detection of Dermatophyte-derived antigen in nails. Scales, scalp specimens, hair or other specimens cannot be used.

- Preparation for specimen collection:** Take a specimen of 1 mg or more according to the guidelines for diagnosis and treatment of cutaneous fungal infection^{1,2}. An inappropriate procedure for specimen collection or an insufficient amount of specimen taken may lead to false negative results or an incorrect judgment. Use clean nail clippers or scissors when taking a specimen and place it in a test tube from the kit.
- Specimen collection must be performed by a professional who is qualified by appropriate education, training and/or experience according to the guidelines for diagnosis and treatment of cutaneous fungal infection. The procedure for specimen collection stated in the guidelines is partially shown below.

Distal and Lateral Subungual Onychomycosis: Remove the area of onycholysis or the tip of the nail, and take a specimen from the deep portion of the nail, as close to the nail bed as possible. If a specimen cannot be taken from the deep portion of the nail, take a specimen from the surface of the skin (actually the nail bed), where onycholysis is present.


Influence of Drugs: The influence of oral antifungals¹ (terbinafine, itraconazole) that are commonly used for the treatment of *tinea unguium* on this assay were evaluated. When the negative control specimen, the positive control specimen and the weakly positive specimen which was prepared by diluting the positive control specimen with the negative control specimen was mixed with the antifungals and subjected to this kit, no influences of these drugs were observed. The concentration of each drug added was approximately 100 times the MIC (minimum inhibitory concentration).

Antifungal	Concentration (µg/mL)	Influence
Terbinafine	0.5	Not observed
Itraconazole	100	Not observed

Test Procedure


Specimen and test device should be at room temperature and the following procedure should be performed at room temperature.

- Prepare the required quantities of test strips, stir rods and extraction buffer. Do not touch the sample pad.

- Place the test strip in the test tube with the sample pad down. Make sure that the sample pad has reached the bottom of the test tube.
- Leave the test strip for at least 5 minutes and determine the result (positive, negative or invalid) by visually checking the presence or absence of coloured bands in the control line zone and the test line zone, within 30 minutes after placing the test strip in the test tube.

Test Interpretation

Precautions for interpretation: If the amount of dermatophyte in the nail specimen is small, the result may become negative. The user should comprehensively evaluate the result of this assay in conjunction with other test results and clinical symptoms. This kit is cross-reactive with other fungi than dermatophyte, such as *Aspergillus* and *Penicillium*. These fungi may be present in the soil or other environments and infect the skin of immune-compromised patients. Professionals should be careful when making a diagnosis.

Test Characteristics

- Results of Clinical Performance Study**⁴: In 222 patients (at 11 centers) suspected of having tinea unguium on visual inspection, a specimen was collected from a foot or hand nail according to the guidelines for diagnosis and treatment of cutaneous fungal infection. The specimen was aliquoted into three portions and subjected to measurement with **PreventID® Dermatophyte**, direct microscopy and PCR (only specimens for which the results of **PreventID® Dermatophyte** and direct microscopy were inconsistent), respectively. Specimen collection, direct microscopy, **PreventID® Dermatophyte** and PCR were performed by different persons under blinded conditions. **Comparison between the results of PreventID® Dermatophyte and PCR incorporating direct microscopy:** Analyses were performed on 222 patients. In 5 patients in whom the results of **PreventID® Dermatophyte** and direct microscopy were inconsistent and PCR could not be performed because the amount of specimen was insufficient, the result of direct microscopy was used.

PCR incorporating direct microscopy				
PreventID® Dermatophyte		Positive	Negative	Total
	Positive	196	5	201
	Negative	6	15	21
	Total	202	20	222

Sensitivity: 97.0%
Specificity: 75.0%
Accuracy: 95.0%

Negative predictive value: 71.4%

Positive predictive value: 97.5%

Comparison between the results of PreventID® Dermatophyte and the dermatologist's final diagnosis (based on the results of direct microscopy, PCR, clinical manifestation, specimen collection site, etc.): Analyses were performed on 217 patients, excluding 5 patients in whom PCR could not be performed because the amount of specimen was insufficient and a final diagnosis could not be made.

Final Diagnosis				
PreventID® Dermatophyte		Dermatophytes	Non-Dermatophytes	Total
	Positive	196	2	198
	Negative	4	15	19
	Total	200	17	217

Sensitivity: 98.0%
Specificity: 88.2%
Accuracy: 97.2%

Negative predictive value: 78.9%

Positive predictive value: 99.0%

- Sensitivity and Accuracy:** When a negative control specimen was tested, this kit provided a negative result. When a weakly positive specimen and a positive control specimen were tested, this kit provided positive results.
- Within-run reproducibility:** When a negative control specimen was tested 4 times, the kit provided a negative result every time. When a weakly positive specimen and a positive control specimen, respectively, were tested 4 times, the kit provided a positive result every time.
- Minimum Detectable Sensitivity:** *Trichophyton rubrum* (NBRC 9185), 0.5 µg dry weight/mL.
- Reference Standard for Calibration:** Dry cells of *Trichophyton rubrum* (NBRC 9185).

- Cross-reactivity:** Autoclaved dry cells of various other fungi than Dermatophyte were added to the extraction buffer at a concentration of 300 µg/mL to evaluate the influence of each fungus on the assay. In addition, colonies of various bacteria grown on agar plates were added to the extraction buffer to evaluate the influence of each bacterium on the assay.

This kit was not reactive with the tested fungi (non dermatophyte) shown below: *Aspergillus nidulans*, *Penicillium citrinum*, *Scopulariopsis brevicaulis*, *Alternaria alternate*, *Pseudallescheria boydii*, *Scedosporium apiospermum*, *Prototheca wickerhamii*, *Schizophyllum commune* (1 nucleus), *Schizophyllum commune*, (2 nucleii), *Absidia corymbifera*, *Basidiobolus ranarum*, *Cunninghamella bertholletiae*, *Mortierella isabellina*, *Mucor circinelloides*, *M. racemosus*, *Rhizomucor pusillus*, *Rhizopus microsporus* var. *rhizopodiformis*, *R. oryzae*, *R. stolonifer* var.*reflexus*, *Syncephalastrum racemosum*, *Zygorhynchus exponents*, *Candida albicans*, *C. dubliniensis*, *C. tolopicalis*, *C. parapsilosis*, *C.guilliermondii*, *C. glabrata*, *C. krusei*, *Geotrichum candidium*, *Trichosporon asahii*, *Cryptococcus neoformans* serotype A, *C. neoformans* serotype B, *C. neoformans* serotype C, *C. neoformans* serotype D, *C. neoformans* serotype AD, *Sporothrix schenckii*, *Fonsecaea pedrosoi*, *Exophiala. jeanselmei*, *Phialophora verrucosa*, *P. richardsiae*, *Rhinocladiella atrovirens*, *Cladophialophora bantiana*, *Malbranchea albolutea*, *M. aurantiaca*, *M. chrysosporioides* *hrysosporioides*, *M. cinnamomea*, *M. dendritica*, *M. filamentosa*, *M. flava*, *M. flocciformis*, *M. fulva*, *M. graminicola*, *M. gypsea*, *M. multicolor*, *M. pulchella*, *Malassezia furfur*, *Gymnoascoideus petalosporus*, *Auxarthron reticulatum*, *Gymnoascus intermedius*, *G. petalosporus*, *G. reessii*, *G. udagawae*, *Emmonsia parva* var. *crescens*, *E. parva* var. *parva*, *Phanerochaete chrysosporium*, *Apinisia queenslandica*, *Arthroderma multifidum*, *Uncinocarpus reesii*, *Chrysosporium carnichaellii*, *C. indicum*, *C. keratinophilum*, *C. pseudomerarium*.

The kit was reactive with the fungi (non dermatophyte) shown below: *Aspergillus flavus*, *A. fumigatus*, *A. niger*, *A. terreus*, *Neosartorya fischeri*, *Paecilomyces lilacinus*, *Penicillium griseofulvum*, *Veronaea botryosa*, *Fusarium solani*, *Exophiala dermatitidis* (M-Y form), *E. dermatitidis* (G form), *E. spinifera*, *Hortaea werneckii*, *Malbranchea circinata*, *M. flavorosea*.

The kit was not reactive with the bacteria shown below: *Escherichia coli*, *Bacillus subtilis*, *Staphylococcus aureus*, *S. faecalis*.

- Reaction with dermatophyte:** Autoclaved dry cells of dermatophyte were added to the extraction buffer at a concentration of 300 µg/mL to evaluate the reactivity of the assay. The kit was reactive with the dermatophyte shown below: *Trichophyton mentagrophytes*, *T. rubrum*, *T. tonsurans*, *T. violaceum*, *T. verrucosum*, *Microsporum gypseum*, *M. canis*, *Epidermophyton floccosum*.

Clinical Significance

PreventID® Dermatophyte, unlike direct microscopy, does not require special skills to determine whether dermatophyte is present or absent³, and this kit, unlike PCR, does not require special equipment. **PreventID® Dermatophyte**, which is easy to use and provides quick results, is an effective assay for the rapid diagnosis of *Tinea unguium*.

Test Limitations

Test results are only reliable if you follow the instructions for use carefully. The test is limited to the detection of dermatophyte in human nails. Use the test only once. Although the test is very accurate, a low incidence of false results can occur. If negative or questionable results are obtained, the test should be repeated on a fresh specimen using a new test device.

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

Références/Naslagwerken/References			
- Ameen M, Lear JT, Madan V, Mohd Mustapa MF, Richardson M. BritishAssociation of Dermatologists’ guidelines for the management of onychomycosis 2014. Br J Dermatol 2014; 171: 937–958. - Guidelines of care for superficial mycotic infections of the skin: onychomycosis. Guidelines/ Outcomes Committee. American Academy of Dermatology. J Am Acad Dermatol 1996; 34: 116–121. - British Association of Dermatologists: Guidelines for treatment of onychomycosis, Roberts DT., Taylor WD., Boyle J., Br. J. Dermatol., 148, 402–410, 2003. - Screening for tinea unguium by Dermatophyte Test Strip, Y. Tsunemi, et al, Br. J. Dermatol., 170, 328–331, 2014. - Clinical study of Dermatophyte Test Strip, an immunochromatographic method, to detect tinea unguium dermatophytes, Y. Tsunemi, M. Hiruma, J Dermatol. 2016 Mar 19. [Epub ahead of print].			
CE	IND		
 Limite de température Temperatuurbegrenzing Temperature limit	 Fabricant Fabrikant Manufacturer	 Contenu suffisant pour <n> tests Inhoud voldoende voor <n> tests Contents sufficient for <n> tests	
 Numéro de lot Chargenummer Batch number	 Numéro de catalogue Catalogusnummer Catalogue number	 Protéger du soleil Uit het zonlicht houden Keep away from sunlight	
 Utilisable jusqu'à Te gebruiken tot Use-by date	 A utiliser avec Te gebruiken met To be used with	 Consulter le mode d'emploi Volg testinstructies Consult instructions for use	
 Ne pas réutiliser Niet herbruikbaar Do not re-use	 Dispositif médical de diagnostic in vitro In-vitro-diagnosemiddel In vitro diagnostic medical device		
 Immundiagnostik Stubenwald-Allee 8a 64625 Bensheim, Germany T: +49 6251 70190-0 info@immundiagnostik.com www.immundiagnostik.com		Distribution/Distributie/Distribution: Preventis GmbH Stubenwald-Allee 8a 64625 Bensheim, Germany T: +49 6251 70711-0 F: +49 6251 70711-299 info@preventis.com www.preventis.com	