REF KST06410G5  
KST06410GP

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### Literatur

1. Shan L, et al.; „Structural basis for gluten intolerance in celiac sprue“; Science; 2002; 297: 2275-9.
2. Comino I, et al.; „Diversity in oat potential immunogenicity: basis for the selection of oat varieties with no toxicity“; Gut; 2011; 60: 915-922.
3. Mariné et al.; „The prevalence of coeliac disease is significantly higher in children compared with adults“; Alimentary Pharmacology Therapeutics; 2011; 33:477-486
4. Silvester J. A., et al.; „Long-term follow-up of individuals with celiac disease: an evaluation of current practice guidelines“; Canadian Journal of Gastroenterology; 2007; 21: 557-564.
5. Walker M. M., et al.; „An update in the diagnosis of coeliac disease“; Histopathology; 2011; 59: 166-179.
6. Comino I, et al.; „Fecal Gluten Peptides Reveal Limitations of Serological Tests and Food Questionnaires for Monitoring Gluten-Free Diet in Celiac Disease Patients“; Am. J. Gastroenterol; 2016;111:1456-1465
7. Sharkey et al.; „Optimization delivery of care in coeliac disease – comparison of the benefits of repeat biopsy and serological follow-up“; Alimentary Pharmacology and Therapeutics
8. Comino I, et al.; „Monitoring of gluten-free diet compliance in celiac patients by assessment of gliadin 33-mer equivalent epitopes in feces“; The American journal of clinical nutrition; 2012; 95: 670-677.
9. Morón B, et al.; „Toward the Assessment of Food Toxicity for Celiac Patients: Characterization of Monoclonal Antibodies to a Main Immunogenic Gluten Peptide“; 2008; PLoS ONE 3:e2294.3.

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REF KST06410G5  
KST06410GP

### Intended Use

**PreventID® GIP Stool** is an immunochromatographic rapid test for detection of gluten immunogenic peptides (GIP) in faecal samples.

### Introduction

Celiac disease is a chronic autoimmune condition characterized by a permanent intolerance of gluten (a mix of proteins present in cereals such as wheat, barley, rye and some kind of oats) [1, 2]. It can arise in genetically predisposed individuals and affects both children and adults. Celiac disease affects an average of 1 in every 133 people, being most frequent in children, where prevalence rises to 1 in 71 [1, 3]. This disease causes the villi in the lining of the small intestine to atrophy, which interferes with the absorption of nutrients such as proteins, fats, carbohydrates, minerals and vitamins. It is a systemic disease that can cause skin lesions, and it also affects joints, the brain and other organs, in addition to the intestine.

Currently, the only effective treatment for celiac disease is to follow a strict lifelong gluten-free diet (GFD). This presents many practical difficulties, especially given that gluten is present in more than 50% of processed non-liquid food products, as an ingredient, an additive or as a consequence of the manufacturing process.

Strict adherence to the GFD is essential to reduce the symptoms of the disease, to avoid nutritional deficiencies and to improve the quality of life of patients. Nevertheless, according to several studies, failure to adhere to the GFD, whether consciously or accidentally, is relatively frequent (32.6 to 55.4%) among celiac patients, especially in adults [4]. Currently there are no direct methods to monitor GFD adherence. The use of antibodies (anti-tissue transglutaminase or anti-deamidated gliadin, whether IgA or IgG) as markers to follow up GFDs is an indirect and generally ineffective method, because it can take months or even more than a year to observe reduced levels of these markers, thus making it impossible to predict recovery on the basis of serological tests. Besides, the levels of these serological markers do not always reflect the state of the intestinal villi [5, 6, 7].

Gluten immunogenic peptides (GIP) are fragments of gluten proteins that are resistant to gastrointestinal digestion and produce immunotoxic reactions in celiac patients. Specifically, they are peptides equivalent to the alpha-gliadin 33-mer, the most immunogenic gluten peptide, that are reactive to the anti-33-mer G12 monoclonal antibody used in this product [8].

The resistance of gluten peptides to gastrointestinal digestion, including the immunogenic 33-mer peptide, ensures that a significant part of the ingested gluten peptides are excreted in feces. Thus, the recovery of measurable quantities of immunogenic fraction of gluten in feces would indicate the consumption of gluten. The presence of gluten immunogenic peptides (GIP) in feces represents a highly useful and direct marker for the accurate short- and long-term control of GFD [9].

### Test Principle

**PreventID® GIP Stool** is an immunochromatographic test that allows the detection of peptides from gastrointestinal degradation of ingested gluten in stool samples. To perform the test, firstly the GIP must be extracted from the faecal samples, being this the most relevant step. Once the sample is extracted, the detection step is based on the reaction of the 33mer-like immunotoxic peptides of gluten in the sample with the coloured conjugates (monoclonal anti-gliadin 33mer antibody/red colored microsphere) previously fixed on the stick [5, 8]. This complex spreads by capillarity through the stick. If the result is positive, a RED line appears in the test zone (T) of the test device. The absence of the RED line indicates a negative result. Whether or not GIP is present, the mixture of the conjugate moves through the stick up to the control region (C) where antibodies have been immobilised and when the test is properly realised a GREEN line (control line) will appear.

### Materials

#### Materials Provided

- Immunochromatographic test devices (with plastic pipette in foil pouch), individually packed [iVYCHECK Sticks and plastic pipette, REF: RS-6405]
- Red-capped tube with extraction solution [Extraction solution, REF: RS-6512]
- Blue-capped tube with dilution solution [Dilution solution, REF: RS-5750]
- Manual

**Material Required but not Provided:** Timer or stop watch, non-powdered disposable gloves, pipette tips, specimen collection container

### Storage and Stability

The product must be stored at a temperature ranging from 2 °C to 30 °C / 35.6 °F to 86 °F during the shelf life of the kit. To obtain optimal test performance, the product must be stored in its original packaging, and used before the expiration date.

NOTE: The foil pouch with the test device should not be opened until its time of use. All components of the kit are fully disposable in ordinary trash or in an appropriate recycling bin.

### Precautions

1. Please read this manual before starting the test. It is recommended to strictly follow the instructions in this manual.
2. To avoid contaminations that interfere with the analysis, the use of non-powdered disposable gloves is recommended. If you do not have disposable gloves, wash your hands thoroughly before the test.
3. Once the **PreventID® GIP Stool** test device has been removed from the foil pouch, it must be used as soon as possible under strict clean conditions.
4. Only for the analysis of human faecal samples.
5. This kit is designed for professional use only.
6. Do not drink any solution (liquid) from the kit. The extraction solution contains alcohol (ethanol).
7. Keep out of reach of children.
8. Do not use any material from the kit after the expiration date.
9. Do not mix components from various kits, or use reagents or different solutions to those supplied.
10. Just touch the test device with gloves or clean hands. In any case, avoid touching the absorbent zone (S), as this may lead to the occurrence of false positives.

### Specimen Preparation and Test Procedure

Stool samples need to undergo an extraction process to make accessible the GIP for further analysis. The complete procedure for the extraction and subsequent analysis of the samples is as follows (also in Fig. 1):

1. Take about 1–2 g or mL (for liquids samples) of faeces. These samples should be collected in a clean and dry container without chemicals, additives or preservatives. Avoid any contact between the sample and urine or between the sample and the water of the toilet. Once the stool samples have been collected, they should be taken to the analytical center within 8 hours of collection. If this is not possible, freeze them at -20 °C. Once the samples arrive at the analytical center: They may be stored at -20 °C for a maximum of 6 months if they arrived frozen or fresh. They can also be analyzed or extracted immediately. If they arrived thawed, samples should be analyzed or, at least, extracted immediately. Sample extracts can be stored at -20 °C for a maximum of two weeks. Faecal samples should be analyzed immediately after being thawed in the laboratory and freeze-thaw cycles should be avoided as they may alter the state of the GIP peptide and cause erroneous results. Once the samples have been extracted and the supernatant transferred to a clean tube, the extracts can be stored for two weeks at -20 °C.


If the fecal sample is frozen, thaw it completely before its analysis.

NOTE: It is recommended that samples are analyzed once completely thawed, avoiding freeze-thaw cycles.

2. Use a clean pipette tip or the spoon attached to the red-capped tube to stir and homogenize the faecal sample. An insufficient homogenization could cause variability in replicates of the same sample.

Take the faecal sample with the spoon attached to the red-capped tube. The correct amount of sample is the one that fills the spoon. Put the spoon with the sample back into the tube and move it to detach the sample from the spoon.

3. Shake the tube containing the sample vigorously intermittently for 10 minutes.
4. Unscrew the small white cap of the red-capped tube and transfer 10 drops of the extracted sample to the blue-capped tube with the dilution solution. Mix by inversion for 15 seconds.
5. Open the foil pouch containing the immunochromatographic test device and use the pipette inside to put 8 drops from the mixture prepared in step 4 onto the S Zone of the test device.
6. The visual interpretation of the results must be done at 10 minutes (a strong positive result will be obtained in samples containing a high amount of GIP). If a RED line is not observed after the first 10 minutes, we strongly recommend to wait an additional 20-minute period (30 minutes in total), since samples containing a low amount of GIP often require more time to show the expected RED line.

### Test Interpretation (Fig. 6)

**Negative:** A single GREEN line (control line) appears in the C zone (Control Zone) of the test device.

**Positive:** In addition to the control line (GREEN), a RED line (Test line) appears in the test zone (T).

**Note:** The intensity of the RED line in the test zone (T) will vary depending on the GIP concentration present in the sample. The presence of any faint line in the Test Zone (T) will indicate positive results no matter how light it is.

**Invalid:** The control line (GREEN) does not appear, whether or not the result line appears (RED). The most common causes for the appearance of an invalid test are: an incorrect procedure, or deterioration of the reagents. In case of invalid results, it is necessary to repeat the experiment with a new test always following a correct procedure. If the problem persists; you must contact the Preventis and stop using the test.

### Sensitivity

The limit of detection of the immunochromatographic test is 0.15 µg GIP/g of faecal sample.

### Specificity

The test is based on the G12 monoclonal antibody which can specifically detect the fraction of wheat prolamins (gliadin), rye (secalin), barley (hordein) that are toxic for celiac patients [1], and also oats (avenin) when in sufficiently high quantities to pose a risk to celiac patients. However, no positive signal is observed when individuals consume rice, corn, buckwheat, soybean, millet, quinoa and amaranth, so the test has no interference with these ingredients free of gluten and thus safe for celiacs.

### Test Limitations

Use the test device only once. Test results are only reliable if you follow the instructions for use carefully. The **PreventID® GIP Stool** is limited to the detection of GIP in human faeces. Although the test is very accurate, a low incidence of false results can occur. If questionable results are obtained, the test should be repeated on a fresh faecal specimen using a new test device.

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

### References

1. Shan L, et al.; „Structural basis for gluten intolerance in celiac sprue“; Science; 2002; 297: 2275-9.
2. Comino I, et al.; „Diversity in oat potential immunogenicity: basis for the selection of oat varieties with no toxicity“; Gut; 2011; 60: 915-922.
3. Mariné et al.; „The prevalence of coeliac disease is significantly higher in children compared with adults“; Alimentary Pharmacology Therapeutics; 2011; 33:477-486
4. Silvester J. A., et al.; „Long-term follow-up of individuals with celiac disease: an evaluation of current practice guidelines“; Canadian Journal of Gastroenterology; 2007; 21: 557-564.
5. Walker M. M., et al.; „An update in the diagnosis of coeliac disease“; Histopathology; 2011; 59: 166-179.
6. Comino I, et al.; „Fecal Gluten Peptides Reveal Limitations of Serological Tests and Food Questionnaires for Monitoring Gluten-Free Diet in Celiac Disease Patients“; Am. J. Gastroenterol; 2016;111:1456-1465
7. Sharkey et al.; „Optimization delivery of care in coeliac disease – comparison of the benefits of repeat biopsy and serological follow-up“; Alimentary Pharmacology and Therapeutics
8. Comino I, et al.; „Monitoring of gluten-free diet compliance in celiac patients by assessment of gliadin 33-mer equivalent epitopes in feces“; The American journal of clinical nutrition; 2012; 95: 670-677.
9. Morón B, et al.; „Toward the Assessment of Food Toxicity for Celiac Patients: Characterization of Monoclonal Antibodies to a Main Immunogenic Gluten Peptide“; 2008; PLoS ONE 3:e22943.

US: all products: Research Use Only. Not for use in diagnostic procedures.

|                            |                                   |
|----------------------------|-----------------------------------|
| Temperature limitation     | Manufacturer                      |
| In vitro diagnostic device | Lot number                        |
| Catalogue number           | Expiry date                       |
| To be used with            | Do not reuse                      |
| Read user instructions     | Contains sufficient for <n> tests |
| Keep away from sunlight    |                                   |

Status: 2019-01-22

Distributed by:

**Preventis GmbH**  
Stubenwald-Allee 8a  
64625 Bensheim, Germany  
Phone: +49 6251 70711-0  
Fax: +49 6251 70711-25  
info@preventis.com  
www.preventis.com