DE Testanleitung

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EN Test Instructions

Intended Use

PreventID® Strep A test is a rapid diagnostic immunoassay for the qualitative detection of Group A Streptococcus antigen. Its analytical sensitivity is 9700 CFU/ml of Strep A cells.

The test is intended as an aid in the rapid diagnosis of a Strep A infection and is for professional use only.

Introduction

Group A Streptococcus (Strep A) is the most common cause of pharyngitis, with the highest morbidity usually found in children. The identification aids in the diagnosis of diseases caused by bacteria belonging to the genus *Streptococcus* and provides epidemiological information on these diseases. Pathogenic streptococci are associated with infections, such as sore throat, impetigo (an infection characterized by small pustules on the skin), urinary tract infections, rheumatic fever, and kidney disease.

Test Principle

The result window of the PreventID® Strep A test cassette is marked with the letters "T" for test line and "C" for control line. No line is visible in the result window before sample application. The control line works as a procedural control. The control line should always appear if the test procedure is performed properly and the test reagents are working properly. Otherwise, the test result is invalid. After sample application, the specimen flows through an absorbing pad and to a conjugate pad containing purified anti-Strep A antibody coupled to purple beads. If the sample contains Strep A antigen, the antibody-coupled to beads will bind to the antigen. As the Strep A antigen-antibody-purple beads are captured on the membrane where purified anti Strep A antibody is spotted, a purple test line (T) will be visible in the result window. If Strep A specific antigen is not detected in the patient sample, no color appears in the test line area (T).

Precautions

- For *in vitro* diagnostic use only.
- Do not eat or smoke while handling specimen.
- Wear protective gloves while handling specimens. Wash hands thoroughly afterwards.
- Avoid splashing or aerosol formation while handling specimen and performing the test.
- Clean up spills thoroughly using an appropriate disinfectant.
- All samples and materials used should be treated as potentially infectious and disposed in a biohazard container. Decontaminate all potentially contaminated objects and surfaces carefully.
- Do not use test if the pouch is torn or if the membrane of the test device is visibly damaged.
- Read the instruction carefully before performing the test.
- Do not mix reagents from different lots.
- For any questions, please contact Preventis GmbH at support@preventis.com.

Storage and Stability

The PreventID® Strep A test kit should be stored at room temperature or refrigerated (4–30°C). The test device is sensitive to humidity as well as to heat. Perform the test immediately after removing the test device from the pouch. Do not use it beyond the expiration date.

Included Materials

- 20 Test cassettes [TEST]
- 20 Tubes
- 1 Instructions for use (IFU)
- 1 Tube rack
- 20 Disposable sample droppers
- 20 Swabs
- 1 Extraction buffer A [BUF A]
- 1 Positive control [CTRL+]
- 1 Extraction buffer B [BUF B]
- 1 Negative control [CTRL-]


Material Required but Not Provided: Timer.

Specimen Collection

- Specimens to be tested directly should be collected by a standard swab method.
- If a sample is to be stored prior to testing, it should be placed in a dry tube, covered, and stored at room temperature or refrigerated.
- All samples should be tested within 5 days after collection.
- If the specimen had been refrigerated, it should be brought to room temperature before testing. Avoid repeated freeze thaw cycles before use.

Specimen Preparation

1. Put the tube into the rack.
2. Put **4 drops extraction buffer A** in the provided tube and then add **4 drops extraction buffer B** to the tube.
3. Place the specimen swab in the tube and swirl it vigorously to mix the re-

agents for about **15 seconds** (Fig. 1). Then incubate the mixture at room temperature for **1 minute** with the swab in the tube.

4. Swirl the swab vigorously for another **15 seconds**, then expunge as much liquid as possible from the swab by pressing and rotating the swab against the wall of the tube. Discard the swab. Mix the contents of the tube by gentle swirling. The mixture is ready for testing.

Test Procedure

1. Remove the test device from the pouch, and place it on a flat, dry surface (Fig. 2).
2. Holding the sample dropper above the test device, squeeze **2 drops** of the mixed specimen into the sample application window. **Note:** If the drops contain many air bubbles, the actual specimen volume may be less than the minimum volume required. So, if there is no red dye migrating to the result window after about 30 seconds, add an additional drop or drops.
3. As the test begins to work, you will see purple colour move across the result window in the center of the test device.
4. Interpret the test results at **10 minutes**.

Caution: The above interpretation time is based on reading the test results at room temperature of 15 to 30°C. All test components and the sample should be brought to this temperature prior testing.

Test Interpretation

- A colour line will appear in the upper section of the result window to show that the test is working properly. This line is the **control line C**.
- The lower section of the result window indicates the test result. If another colour line appears in the lower section of the result window, this line is the **test line T**.

Positive: The presence of **two colour lines** (T and C) within the result window (regardless of which line appears first) indicates a positive result. **Note:** Generally, the higher the analyte level in the specimen, the stronger the colour of the test line will be. When the specimen analyte level is close to but still within the sensitivity limit of the test, the colour of the test line (T) will be very faint.

Negative: The presence of only one colour line (C) within the result window indicates a negative result.

Invalid: If after performing the test, no colour lines are visible within the result window, the result is considered invalid. The test result is also considered invalid if only the test line (T) appears. Some causes of invalid results are: not following the directions correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen should be re-tested using a new test kit.

Note: A positive result will not change once it has been established at 10 minutes. However, in order to prevent any incorrect results, the test result should not be interpreted later than after 10 minutes.

Quality Control

External controls may be tested when opening a new test kit.

1. Add 2 to 3 drops of positive control or negative control solution into the sample application window of the test device. Interpret test results at 10 minutes. Do not interpret test results later than after 10 minutes.
2. Please see the previous section "Test Interpretation" for interpreting the test results.

Limitations


Although PreventID® Strep A is very accurate in detecting Strep A, a low incidence of false results can occur. Other clinically available tests are required if questionable results are obtained. PreventID® Strep A is a qualitative assay. The amount of Strep A present in the specimen cannot be determined by this test. The results merely distinguish positive from negative samples. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

Expected Values

Approximately 15% of pharyngitis in children ages 3 months to 5 years is caused by Group A beta-hemolytic *Streptococcus*.⁶ In school-aged children and adults, the incidence of Strep throat infection is about 40%.⁷ This disease usually occurs in the winter and early spring in temperate climates.

Test Performance Sensitivity and Specificity

A total of 310 throat swabs were collected from patients exhibiting symp-

toms of pharyngitis. Each swab was rolled onto a sheep blood agar plate, and then tested by PreventID® Strep A. The plates were further streaked for isolation, and then incubated at 37°C with 5–10% CO₂ and a Bacitracin disk for 18–24 hours. The negative culture plates were incubated for an additional 18–24 hours. Possible GAS colonies were subcultured and confirmed with a commercially available latex agglutination grouping kit. Of the 310 total specimens, 160 were confirmed to be negative and 150 were confirmed to be positive by culture. During this study, two Strep F specimens yielded positive results with the test. One of these specimens was re-cultured, then re-tested and yielded a negative result. Three additional different Strep F strains were cultured and tested for cross-reactivity and also yielded negative results.

	Method		Total Results
	+	-	
PreventID® Strep A	+	143	153
	-	7	157
Total Results	150	160	310
Accuracy	95.3%	93.7%	

Cross-Reactivity

The following organisms were tested at 1.0 × 10⁷ organisms per test and were all found to be negative when tested with the Strep A test device. No mucoid-producing strains were tested: *Branhamella catarrhalis* (ATCC 43628), *Candida albicans* (ATCC 14053), *Corynebacterium diphtheriae* (ATCC 9015), *Escherichia coli* (ATCC 25922), *Haemophilus influenzae* (ATCC 35056), *Klebsiella pneumoniae* (ATCC 13883), *Neisseria gonorrhoea* (ATCC 9826), *Neisseria meningitidis*, serogroup B (ATCC 13090), *Neisseria sp.* (ATCC 43831), *Pseudomonas aeruginosa* (ATCC 27853), *Serratia marcescens* (ATCC 8100), *Staphylococcus aureus* (ATCC 29213 & 25923), *Staphylococcus epidermidis* (ATCC 12228), *Streptococcus Group B* (ATCC 12386), *Streptococcus Group C* (ATCC 12388), *Streptococcus Group D* (ATCC 12389), *Streptococcus Group F* (ATCC 12393), *Streptococcus Group G* (ATCC 12394), *Streptococcus pneumoniae* (ATCC 9163, 6306 & 10015).

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Gebrauchsanweisung beachten
Consult instructions for use

Verwendbar bis
Use-by date

In-vitro-Diagnostikum
In vitro diagnostic medical device

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Katalognummer
Catalogue number

EU REP
EU-Bevollmächtigte
Authorized representative in EU

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Ameritech Diagnostic Reagent (Jiaxing) Co., Ltd.
K4-2 Science Technology Garden, Economic Development Zone
314500 Tongxiang, Zhejiang, P. R. China

Vertrieb / Distribution:

Preventis GmbH
Stubenwald-Allee 8a, 64625 Bensheim, Germany
T: +49 6251 70711-0 • F: +49 6251 70711-299
info@preventis.com • www.preventis.com

AR Experts BV
Boeingavenue 209
1119 PD Schiphol-Rijk
The Netherlands
www.ar-experts.eu