DE Testanleitung

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EN Test Instructions

Intended Use

PreventID® Strep B is a rapid test for the detection of group B streptococci from swabs or culture. Its analytical sensitivity is 4500 CFU/mL of Strep B cells. The test is for professional use only.

Introduction

The identification aids in the diagnosis of diseases caused by bacteria belonging to the genus *Streptococcus* and provides epidemiological information on these diseases. Pathogenic streptococci are associated with infections, such as sore throat, impetigo (an infection characterized by small pustules on the skin), urinary tract infections, rheumatic fever, and kidney disease.

Test Principle

The result window of the PreventID® Strep B test cassette is marked with the letters "T" for test line and "C" for control line. No line is visible in the result window before sample application. The control line works as a procedural control. The control line should always appear if the test procedure is performed properly and the test reagents are working properly. Otherwise, the test result is invalid. After sample application, the specimen flows through an absorbing pad and to a conjugate pad containing purified anti-Strep B antibody coupled to purple beads. If the sample contains Strep B antigen, the antibody-coupled to beads will bind to the antigen. As the Strep B antigen-antibody-purple beads are captured on the membrane where purified anti Strep B antibody is spotted, a purple test line (T) will be visible in the result window. If Strep B specific antigen is not detected in the patient sample, no color appears in the test line area (T). The test's accuracy does not depend on the organism's viability. Instead, group B strep antigen is extracted directly from the swab and identified using antibodies specific for the group B strep carbohydrate.

Precautions

- For *in vitro* diagnostic use only.
- Do not eat or smoke while handling specimen.
- Wear protective gloves while handling specimens. Wash hands thoroughly afterwards.
- Avoid splashing or aerosol formation while handling specimen and performing the test.
- Clean up spills thoroughly using an appropriate disinfectant.
- All samples and materials used should be treated as potentially infectious and disposed in a biohazard container. Decontaminate all potentially contaminated objects and surfaces carefully.
- Do not use test if the pouch is torn or if the membrane of the test device is visibly damaged.
- Biotin can interfere and cause incorrect test results for subjects who have elevated levels of biotin in their specimen.
- Read the instruction carefully before performing the test.
- Do not mix reagents from different lots.
- For any questions, please contact Preventis GmbH at support@preventis.com.

Storage and Stability

The PreventID® Strep B test kit should be stored at room temperature (4–30°C). The test device is sensitive to humidity as well as to heat. Perform the test immediately after removing the test device from the pouch. Do not use it beyond the expiration date.

Included Materials

- 20 Test cassettes [TEST]
- 1 Instructions for use (IFU)
- 20 Disposable sample droppers
- 2 Bottles with extraction buffer [BUF]
- 20 Swabs
- 20 Tubes
- 1 Tube rack
- 1 Positive control [CTRL+]
- 1 Negative control [CTRL-]

Material Required but Not Provided: Timer.

Specimen Collection

- Swab the lower vagina (vaginal introitus), followed by the rectum (i.e., insert swab through the anal sphincter) using the same swab or two different swabs.
- Do not use swabs with cotton or calcium alginate tips or wooden shafts.
- Do not use swabs impregnated with charcoal or transport media containing agar gelatin.
- If a sample needs to be stored prior to testing, it should be placed in a dry tube, covered, and refrigerated.
- All samples should be tested within 5 days after collection.
- If specimen had been refrigerated, it should be brought to room temperature before testing. Avoid repeated freeze thaw cycles before use.

Specimen Preparation

1. Put the tube into the rack.
2. Pipet **12 drops extraction buffer** into the tube.
3. Place the specimen swab in the tube and swirl it vigorously to mix the reagents for about **15 seconds** (Fig. 1). Then incubate the mixture at room temperature for **10 minutes** with the swab in the tube.

4. Swirl the swab vigorously for **15 seconds**, then expunge as much liquid as possible from the swab by pressing and rotating the fiber portion against the wall of the tube. Discard the swab. Mix the contents of the tube by gentle swirling. The mixture is ready for testing.

Test Procedure

1. Remove the test device from the pouch, and place it on a flat, dry surface (Fig. 2).
2. Holding the sample dropper above the test device, squeeze **3 drops of the mixed specimen** into the sample application window (Fig. 2). Wait until each drop is absorbed, before adding additional drops. **Note:** If the drops contain many air bubbles, the actual specimen volume may be less than the minimum volume required. So, if there is no red dye migrating to the result window in about 30 seconds, add an additional drop or drops.
3. As the test begins to work, you will see a purple line move across the result window in the center of the test device.
4. If the test does not begin to work (no purple line moving across the result window) after the application of the specimen, add 1 to 2 drops of buffer directly from the buffer bottle.
5. Interpret the test results at **10 minutes**. Do not interpret test later than after 10 minutes.

Caution: The above interpretation time is based on reading the test results at room temperature of 15 to 30°C. All test components and the sample should be brought to this temperature prior testing.

Test Interpretation

- A colour line will appear at the left section of the result window to show that the test is working properly. This line is the **control line (C)**.
- The right section of the result window indicates the test result. If another colour line appears at the right section of the result window, this line is the **test line (T)**.

Positive: The presence of **two colour lines** (T and C) within the result window (regardless of which line appears first) indicates a positive result. **Note:** Generally, the higher the analyte level in the specimen, the stronger the colour of the test line (T) will be. When the specimen analyte level is close to but still within the sensitivity limit of the test, the colour of the test line (T) will be very faint.

Negative: The presence of only one colour line (C) within the result window indicates a negative result.

Invalid: If after performing the test, no colour lines are visible within the result window, the result is considered invalid. The test result is also considered invalid if only the test line (T) appears. Some causes of invalid results are: not following the directions correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen should be re-tested using a new test kit.

Note: A positive result will not change once it has been established at 10 minutes. However, in order to prevent any incorrect results, the test result should not be interpreted later than after 10 minutes.



Quality Control

External controls may be tested when opening a new test kit. External controls may be performed with the positive and negative control provided in the test kit.

1. Add 2 to 3 drops of positive control or negative control solution into the sample application window of the test device. Interpret test results at 10 minutes. Do not interpret test results later than after 10 minutes.
2. Please see the previous section "Test Interpretation" for interpreting the test results.

Limitations

Although the test is very accurate in detecting Strep B, a low incidence of false results can occur. Other clinically available tests are required if questionable results are obtained. As is true with any diagnostic procedure, the physician should evaluate data obtained by the use of this kit in light of other clinical information, including culture, if results are inconsistent with clinical symptoms. The PreventID® Strep B is a qualitative assay. The amount of Strep B present in the specimen cannot be estimated by this as-

say. The test results distinguish positive from negative samples. A positive result indicates the sample contains Strep B above the cut-off concentration. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

Test Performance

Relative Sensitivity and Specificity

When compared with ELISA test, PreventID® Strep B test has a **relative sensitivity** of 97.4% and a **relative specificity** of 98%.

Literatur / References

1. Brady K. et al. (1989). Reliability of a Rapid Latex Fixation Test for Detecting Group B Streptococci in the Genital Tract of Parturients at Term. *Obstet Gynecol*, 73(4), 678–681.
2. Stillier R. J. et al. (1989). Rapid Detection of Vaginal Colonization With Group B Streptococci by Means of Latex Agglutination. *Am J Obstet Gynecol*, 160(3), 566–568.

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Gebrauchsanweisung beachten
Consult instructions for use

Verwendbar bis
Use-by date

In-vitro-Diagnostikum
In vitro diagnostic medical device

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Katalognummer
Catalogue number

EU REP
EU-Bevollmächtigte
Authorized representative in EU

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Ameritech Diagnostic Reagent (Jiaxing) Co., Ltd.
K4-2 Science Technology Garden, Economic Development Zone
314500 Tongxiang, Zhejiang, P. R. China

Vertrieb / Distribution:

Preventis GmbH
Stubenwald-Allee 8a, 64625 Bensheim, Germany
T: +49 6251 70711-0 • F: +49 6251 70711-299
info@preventis.com • www.preventis.com

AR Experts BV
Boeingavenue 209
1119 PD Schiphol-Rijk
The Netherlands
www.ar-experts.eu