de

PREVENTIS

PREVENTIS *SmarTest*[®]

FIT Home

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PST40111EP

CE 0483

Immunological rapid test for self-testing for the quantitative detection of occult blood in stool as an aid in the check-up for colon cancer

Test instructions

REF PST40111EP

played in the app. Carry out the camera test at the same location where you want to evaluate the test later on to include the local lighting conditions into the check. Please note: If the camera test repeatedly fails, check the camera lens for scratches and contaminants and remove them. Check as well whether the camera test card is eventually damaged and take care for an even illumination. If the test continues to fail, please contact the technical service of Preventis.

- Specify a PIN. Follow the instructions on your smartphone. First enter a PIN and confirm it in the next step by repeating the assigned PIN. Caution: The PIN is used to protect your data and cannot be restored.
- Specify whether you want to use the reading function of the test instructions. You can change this setting at any time in the app settings.
- The **SmarTest FIT**[®] app is now set up and ready to use right away.

Functions of the SmarTest FIT[®] app

In the main menu of the **SmarTest FIT**[®] app the following functions are available: "Prepare test" and "Results". Via „Prepare test“ you reach the menu for the guided test procedure. Images, texts and animations explain the test procedure.

Hence you can follow the test well, even if you are not yet familiar with the test or have not performed it for some time. "Results" gives you an overview about all results measured with your smartphone. The „gear“ symbol takes you to the settings.

Before performing the test, please read the following instructions.

Sample Preparation

Note:

- To avoid false positive results, the sample should not be taken during or up to 3 days after menstruation and not with severe bleeding gums, bleeding haemorrhoids or blood in urine. Also samples with clearly visible blood in the stool or reddish coloured stool extracts are not suitable for analysis with **SmarTest FIT**[®] Home and must not be used.
- No dietary restrictions are necessary before the test.
- Alcohol, aspirin and some medications sometimes cause gastrointestinal irritation or bleeding. Ask your doctor prior to the test to prescribe a temporary dose reduction if necessary.

Sample Collection (Fig. 1)

- Gain the stool sample using the paper stool catcher according to the instructions. The stool sample must not come in contact with water or urine.

Test Characteristics

Analytical Sensitivity: The test has a limit of detection (LOD) of 5.4 µg haemoglobin/g stool. The measuring range is from 5.4 to 70 µg/g.

Analytical Specificity: The test shows no cross-reactivity with haemoglobin of animal origin.

The deviations in the measurement of different batches have coefficients of variation of ≤10%. The precision of the test (intra-assay accuracy under constant conditions) shows a coefficient of variation of ≤22.2%. For reproducibility (inter-assay accuracy with different conditions for person, date, smartphone, batch), the coefficient of variation is ≤28.1%.

Cross-reactivity: No interference was found testing the following substances/concentrations. Testing was done by adding 1/100 of the depicted concentration to an extracted stool sample: soy bean oil 30%, mucin 5%, Diclofenac 0.1%, and cyclamat-saccharin-mix 1.3%. Evaluation was done using 10 replicates for each tested concentration/substance in comparison with the control group.

Clinical Sensitivity and Specificity: With a cut-off at 9.6 µg haemoglobin/g stool **SmarTest FIT**[®] Home shows a clinical specificity of 93.0% and a clinical

sensitivity of 31.5% for advanced neoplasia (colorectal carcinoma and large adenomas). With this cut-off the expected positive rate in an organised colon cancer screening (= screened representative selection of collected stool samples) is 10.4% [Gies et al., 2018].

Materials

Materials provided

- 1 test device, individually packed **TEST** (with desiccant, not required for testing)
- 1 sample collection tube, including 1.5 ml buffer solution and sample collection stick **TUBE**
- 1 paper stool catcher
- 1 camera test card **CARD**
- test instructions

Additional materials required: disposable gloves

Note: In order to obtain the rapid test results you will need the **SmarTest FIT**[®] app, which can be downloaded from the Apple App Store or Google Play Store.

Storage and Stability

SmarTest FIT[®] Home test system should be stored between 4–30 °C. Do not use the test components beyond the expiry date. Protect the test from heat and direct sunlight.

General Notes

Please allow around 20 minutes to run the test without interruption. Ensure that your smartphone has sufficient power (at least 20% battery level). Using **SmarTest FIT**[®] app requires an active Internet connection at regular intervals, at least every 48 hours, to download the latest calibration data. While running the test, do not switch off your smartphone or switch it to silent as the timer for the test will emit a beep.

Precautions

- For *in vitro* diagnostic use only. Only for application outside the body. For single use only.
- Treat the sample as potentially infectious: do not eat or smoke while handling specimen. Wear protective gloves if available and wash hands thoroughly.
- Do not use the test after the expiry date has passed.
- Do not use the test if the packaging or test components are damaged.
- Do not drink the content of the sample tube and avoid direct contact with the liquid. If some of the liquid comes into contact with your skin, thoroughly wash the affected area immediately.
- The test is designed to be used at room temperature and moderate humidity (<60 % relative humidity). Do not use the test with high humidity or temperature.

- Only remove the test device from the aluminium pouch immediately before running the test.
- Keep out of the reach of children.

- The product quality is constantly monitored. If you have any questions please contact Preventis GmbH (info@preventis.com).

Installing the SmarTest FIT[®] app

Please note: The installation of the **SmarTest FIT**[®] app is mandatory for the evaluation of the rapid test!

The **SmarTest FIT**[®] app is available for:

- iOS, from the Apple App Store
- Android, from Google Play Store

The already supported smartphones (system requirements) can be found here: [fit.preventis.com](https://www.preventis.com).

If your smartphone is not yet listed, you can use the camera test card to qualify it for the test. **Note:** not all smartphones meet the requirements of the app.

Setting up the SmarTest FIT[®] app

Before initial use, the **SmarTest FIT**[®] app must be set up. To do this, carry out the following:

- Start the **SmarTest FIT**[®] app on your smartphone.
- Read the General Terms and Conditions carefully and confirm that you agree to them by clicking on "continue".
- Carry out the camera test by using the camera test card. To do this, follow the instructions dis-

- Keep the sample collection tube – subsequently called TUBE – upright when opening it so that no buffer solution leaks out. The correct amount of buffer solution is important for the accuracy of the result. Carefully handle the tip of the TUBE so that it does not break off.

Hold the yellow sample collection stick at the top and slowly pull it out of the TUBE by turning it counterclockwise. Stick the tip of the sample collection stick in **at least 3 different locations in the stool sample in one go**.

Attention: All grooves of the sample collection stick must be filled with stool. Gently wipe off excess stool from the tip of the sample collection stick on the paper stool catcher so that only the amount of stool sticking in the grooves gets into the TUBE.

- Put the sample collection stick with the stool sample **only once** back into the TUBE filled with buffer solution and do not pull it out.

Attention: The buffer solution in the TUBE changes its colour hardly. Repeated transfer of stool has to be avoided because your sample otherwise can not be measured or false results will be obtained.

- Close the TUBE by turning the sample collection stick clockwise. Shake vigorously the closed TUBE until the stool sample has detached from the sample collection stick.

- The paper stool catcher can be flushed down the toilet.
- After sampling preparation perform the test as soon as possible. If necessary, the tube can be stored at room temperature for **up to 3 days**.


Fig. 1: Sample collection

Test Procedure and Test Interpretation

Have the stool sample collection tube filled with stool extract plus the test device (packed) ready. Start the **SmarTest FIT**[®] app on your smartphone.

Test procedure: „Prepare test“

Here you get an illustrated and animated test instruction directly on your smartphone. Start the guided test procedure by pressing the „Prepare test“ button in the main menu of the **SmarTest FIT**[®] app. Before the actual test procedure begins, you will be guided through some preparatory steps:

- Please take sufficient time to perform the test. The execution time is approximately 20 minutes.
- Make sure that all components of the test kit are at room temperature. Even if the **SmarTest FIT**[®] Home tests are stored at room temperature, it is possible that the test components have a different temperature than the surrounding room, for example due to transport.
- Check your smartphone's Internet connection. **SmarTest FIT**[®] requires an active Internet connection within 48 hours since the last test to download the latest calibration data. In order to inform you about the incubation time by means of an alarm, the ringtone of the smartphone must be set to loud.

- In the app you have the possibility to call up further information from the test instructions by clicking on the i-symbol during the test execution.

Remove the test device from its packaging and place it on a flat, dry surface.

Take the sample collection tube into your hand. Now shake the sample collection tube vigorously, hold it with the tip upwards, wrap a soft piece of paper around the tip and carefully break off the tip.

Hold the sample collection tube approx. 1 cm above the sample application window and place **3 hanging drops of the sample** into the sample application window of the test device by gently pressing the sample collection tube (Fig. 2). It is important that the drops fall freely to ensure the correct drop volume.


Fig. 2: Test device

After applying the sample, immediately press the „Start test“ button in the **SmarTest FIT**[®] app. A timer checking the incubation time starts automatically. You will hear a beep after the 15-minute incubation period has elapsed. You have now 5 minutes to scan the result. If this time is exceeded, no guarantee can be given for the correctness of the result! In that case the measurement has to be repeated with a new test set.

Test interpretation

First check to see whether there is a red control line in the result window of the test device. The control line in the result window should be located underneath the letter "C" on the test device. If no control line is visible, the test is invalid and cannot be analysed. The position of the test line is at "T" (Fig. 3).


Fig. 3: Test interpretation

To scan the rapid test place the test device on a bright and flat surface. Make sure the light conditions are similar to those in the camera test. Avoid shadows, strong light from the side and direct sunlight.

Then select in the **SmarTest FIT**[®] app "Scan Test". You will see a template of the test device on your screen. Align the test device in front of you with the template on the screen. While doing this, keep your smartphone horizontally parallel to the test device and avoid tilting it. Hold this position briefly until the **SmarTest FIT**[®] app automatically triggers the camera and goes to the analysis screen.

After the calculation the result will automatically be displayed directly on the screen of your smartphone. Press "Continue" to proceed to the subsequent screen for result explanation and recommendations on follow-up actions.

Your measurement is finished now. Use the button "Continue" to return to the main menu.

Interpretation of the Result

Positive: haemoglobin concentrations ≥9.6 µg/g stool (more or equal) are interpreted as positive, meaning human blood was detected in a pathological amount. Blood in the stool can have many causes and is only in some cases due to bleeding intestinal polyps or tumours. Other common causes include menstrual bleeding, heavy bleeding of the gums, haemorrhoids, stomach irritations or blood in the urine. The evidence of blood provided here is therefore not directly equivalent to a medical diagnosis of colon cancer. **SmarTest FIT**[®] Home is just one test in a necessary series of tests to enable your doctor to make a complete diagnosis. We therefore recommend that you **visit your doctor as soon as possible** in order to better classify the results of the rapid test.

Negative: haemoglobin concentrations <9.6 µg/g stool (less than) are interpreted as negative, meaning no human blood was detected in a pathological amount. Even if your stool sample today did not show an elevated level of occult blood, bleeding cannot completely be excluded. Occult blood is not distributed in the stool evenly. Colon cancer and colon polyps, especially in the early stages, do not bleed at all times.

To increase your safety, repeat the measurement regularly – preferably once a year – or talk to your doctor about further preventive measures.

Test Limitations

The **SmarTest FIT**[®] Home test system detects occult blood with a high degree of reliability; however, in certain cases it may produce incorrect results. In case of questionable results, other clinical tests are indicated. As with all diagnostic tests, a definite clinical diagnosis should not be based on a single test result, but should be made by a qualified professional after evaluation of all clinical and laboratory findings.

Disposal

After use please dispose all parts of **SmarTest FIT**[®] Home test set with the residual waste. You can flush down the paper stool catcher in the toilet.

Literature:
Gies A, Cuk K, Schrotz-King P, Brenner H (2018) Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening. *Gastroenterology* 154(1):93-104. doi: 10.1053/j.gastro. 2017.09.018. Epub 2017 Sep 25.

 Temperature limitation	 Manufacturer
 In vitro diagnostic device: only to be used outside of the body	 Lot number
 Catalogue number	 Expiry date
 To be used with	 Do not reuse
 Read user instruction	 Contains sufficient for <n> tests
 Keep away from sunlight	

2025-03-26