PREVENTIS *SmarTest*[®]

Ferritin Home

Testanleitung
Instructions for use
Mode d'emploi
Istruzioni per l'uso
Kullanma Kılavuzu

PREVENTIS

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Ferritin Home

Intended purpose

Preventis SmarTest® Ferritin Home is a rapid test for the quantitative immunological determination of the iron storage protein Ferritin in capillary blood. This value aids in diagnosis of absolute or relative (functional) iron deficiency or overload. The manual test system consists of an antibody-based strip test, a so-called lateral flow test, including a smartphone app developed for test evaluation. The test is an *in vitro* diagnostic and suitable for self-testing by laypersons 18 years of age and older. **The package contains materials to perform one, single-use test.**

General information

- **Preventis SmarTest® Ferritin Home** is intended for *in vitro* diagnostic use only, which means for use outside of the body only.
- Each test pack contains materials for **one, single-use test**.
- Test components from different production units must not be mixed.
- All test components may only be used for the designated purpose and in combination with the components included in the kit.
- The prepared sample may only be used in combination with the associated test components.

- Performance of the test does not require any special expertise and can be carried out at home. Please allow at least 20 minutes to perform the test undisturbed.
- The test must be performed at room temperature and a humidity of 20–80%.
- To correctly classify your test result, the specification of your biological sex is necessary.
- In the presence of infection or inflammatory events, the ferritin value may be elevated. Please do not perform the test until symptoms are no longer present.
- Read test instructions carefully before use. Set up the app, take the sample, perform the test and evaluate it exactly according to the test instructions.
- Make sure your smartphone is sufficiently charged (battery level at least 20%).
- The **SmarTest Ferritin** app requires an active internet connection at regular intervals, but at least every 48 hours, to load the latest calibration data.
- Do not turn off the phone and keep the ringtone on during the application, as the timer uses a signal tone.
- Perform the device qualification and camera test in the same location that the rapid test will be performed in order to accurately simulate lighting conditions. Be sure that the lighting is even. Avoid shadows, light from the sides and direct sunlight.
- Ensure that the camera test card and test cassette are placed on a flat, dry and bright surface, resting horizontally at all times.
- The chase buffer must be added to the test cassette 30 seconds after sample application. You have a maximum of 1 minute to do this. You will only receive a valid result if the chase buffer is

applied within this time period.

- After the 15-minute timer has expired, you have 2 minutes to scan the test cassette with your smartphone. You will only receive a valid result if the test is evaluated within this time period.
- For the security of your data, it will only be transmitted in encrypted form. For exact details please refer to the corresponding Data Privacy Policy from **Preventis GmbH**, which can be found at www.preventis.com.
- Please use only the latest version of your original operating system provided by the manufacturer of your smartphone. A list of supported operating systems can be found at link.preventis.com/ferritin-home.
- Ensure that the most recent version of the app is installed.

Safety notes

- Do not drink the contents of the chase buffer, avoid direct contact with the liquid. Immediately wash contaminated areas thoroughly after skin contact.
- Keep all components out of reach of children.
- If you have problems setting up the app, collecting the specimen, performing or evaluating the test, contact info@preventis.com or have a third person, e.g. a family member, help you.
- Treat the specimen as potentially infectious: do not eat, drink, or smoke during specimen handling. It is recommended to always wash/disinfect hands after testing. **In case of serious incidents, please contact Preventis GmbH (info@preventis.com) and the appropriate authorities immediately.**

Storage, stability and disposal

- Store the test at room temperature or in the refrigerator (+4 to +30°C), do not freeze.
- The test is sensitive to direct sunlight and high temperatures.
- The test must be used immediately after opening the test cassette pouch (within 5 minutes).
- The test tolerates a humidity of 20–80%.
- Do not use the test if the packaging or the supplied materials are damaged.
- Do not use the test after the expiration date.
- After use, dispose of all test components, including the camera test card and test instructions, in the household waste.

Medical background^{1,2,3,4,5}

Iron is needed for the formation of the red blood pigment hemoglobin, however, in its free form it is toxic to the body. The free iron molecules are bound by the protein ferritin and are thereby stored. Ferritin is present in every cell of the body and is found predominantly in the liver, spleen and bone marrow. The ferritin content in the blood serves as a measurement for the size of the body iron stores and can indicate an oversupply or deficiency. Ferritin is the best marker for detecting iron deficiency at an early stage.

Regular monitoring of the ferritin status is especially recommended for people who are susceptible to iron deficiency or excess. This includes vegans/vegetarians, people with celiac disease, Crohn's

disease, iron storage disease or liver disease, women with heavy menstrual periods, pregnant women, and regular blood donors or athletes. Testing of iron status is also recommended in case of symptoms such as dizziness, fatigue, reduced ability to perform and concentrate, headache, pallor or weight loss. Ferritin levels may be elevated in cases of acute or chronic inflammation in the body or chronic liver disease. During pregnancy and lactation, the iron requirement is increased.

Test principle

Preventis SmarTest® Ferritin Home is an antibody-based strip test for the detection of ferritin in capillary blood. The sample is collected using a specific sampling system and then applied to the test cassette. Now the timer is started. After 30 seconds a signal sounds and 1 drop of the chase buffer is placed on the sample application window (S) of the rapid test.

Ferritin present on the test strip reacts with gold-labeled anti-ferritin antibodies and – if performed properly – becomes visible as a reddish test band after 15 minutes at the latest. The appearance of the control band (C) also indicates whether the test has functioned correctly from a technical point of view (see **Evaluation**, page 36). Using the **SmarTest Ferritin** app, the quantitative result is then evaluated based on the color intensity of the control and test bands. The collection of the sample before the measurement, as well as the test procedure and evaluation are described step by step in the chapters **Sample extraction** (see page 31), **Test procedure** (see page 33) and **Evaluation** (see page 36).


Included materials

1 test cassette [TEST], individually wrapped with desiccant bag*

*To prevent the test strip from absorbing moisture, a desiccant bag is included in the packaging. Dispose of the desiccant bag after opening.

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[CHASE]

Additional materials needed

A smartphone with the **SmarTest Ferritin** app and an internet connection.

Setting up the SmarTest Ferritin app

The **SmarTest Ferritin** app guides you step-by-step through the sampling and testing process and displays your test result at the end.

Download and installation

Attention: The successful installation (app icon appears on your smartphone) of the **SmarTest Ferritin** app is mandatory for the evaluation of the quick test! The **SmarTest Ferritin** app is available for **iOS**, available in the App Store and **Android**, available in the Google Play Store.

You can find the already supported smartphone models and operating systems at: link.preventis.com/ferritin-home. If your smartphone is not yet listed, you can qualify it for the test using the camera test card (CARD) (see **Device qualification**, page 29).

Download from
App Store here

Download from
Google Play Store

Please note: Unfortunately, not all smartphones meet the requirements of the app. If you have any questions, please feel free to contact us at info@preventis.com.

If you do not receive the device qualification prompt, then go to the next step **Camera test**, page 30. Your smartphone model has then already been approved by us.

mySmarTest

This is an optionally usable portal where your test results can be saved. **It provides you with an overview of your measurements to date.** In addition, you can use **mySmarTest** to share test results with healthcare professionals via the „**Share results**“ menu item if you wish (see **Share results**, page 28). You can set up your account directly during app setup or via „**mySmarTest**“ on the main screen.

Account information

Your login data consist of your e-mail address and a password that you assign yourself. **These login data are stored independently of your test results.** If you forget your password, your account can be recovered using your email address.

Share results (optional)

Open the **SmarTest Ferritin** app and press „**mySmarTest**“ on the main screen. You can then share your test results with your treating healthcare professional via „**Share Results**“ if you wish. You can adjust the access rights at any time under „**Manage shared results**“.

Device qualification

Attention: Perform the camera test in the same location that the rapid test will be performed later in order to accurately simulate the lighting conditions. Be sure that the lighting is even. Avoid shadows, light from the sides and direct sunlight.


Step 1

Turn the camera test card (CARD) in such a way that the green stripe is on the right side.

Note: If the app does not prompt you to perform the device qualification, go directly to the **Camera test** chapter, page 30.

Step 2 (Fig. 2, page 29)

Align the template on your display with the six squares on the camera test card. If the alignment is correct, the template will turn green. Make sure to hold your smartphone parallel to the camera test card until you receive a success notification.

Note: Please make sure to hold your smartphone horizontally to the camera test card lying in front of you. Vary the distance of your smartphone to the camera test card until the templates on your screen match the outlines of the squares.

If the device qualification fails repeatedly, check the camera lens for scratches and dirt

and remove them if possible. Also check the camera test card for any damage.


Attention: If the device qualification fails repeatedly, your smartphone may not be suitable for evaluation of our rapid tests. For further questions please contact us at info@preventis.com.

Camera test

Attention: Perform the camera test in the same location that the rapid test will be performed later in order to accurately simulate the lighting conditions. Be sure that the lighting is even. Avoid shadows, light from the sides and direct sunlight.

Step 1 (Fig. 1, page 30)

Place the camera test card in front of you in such a way that the green stripe is facing up and you can read the numbers. Align the template on your display with the image of test cassette ① until you receive a success message.

Step 2

Repeat this step with the image of test cassette ②.

Note: Make sure to hold your smartphone parallel to the camera test card. If the camera test fails repeatedly, check the camera lens for scratches and dirt and remove them if possible. Also check the camera test card for any damage.

Please make sure to hold your smartphone horizontally to the camera test card lying in front of you. Vary the distance of your smartphone to the camera test card until the template on your screen matches the outline of the depicted test cassette. If the test continues to fail, your smartphone may not be suitable for evaluation of our rapid tests. For further questions please contact us at info@preventis.com.

Sample extraction

Note: Please allow for approximately 20 minutes to complete the test. Make sure that the test materials are at room temperature. Get the sealed test cassette and all other test components in place. Make sure that the tip of the capillary does not come into contact with anything.

Attention: Leave the test cassette in its packaging. Check the internet connection of your smart-

phone and turn on the sound. Open the **SmarTest Ferritin** app and select „**Prepare test**“. Then start the animated video instructions for sample extraction, test procedure and evaluation via „**To test**“.

Please perform all steps one after another without pauses. Do not store the sample in the capillary, but always immediately proceed with the next steps. After sample application, add the chase buffer to the test cassette immediately after the 30 seconds have elapsed and the signal in the **SmarTest Ferritin** app has sounded.

Step 1

The lateral areas of the fingertip of the middle or ring finger are particularly suitable for collecting capillary blood. Rub the selected finger until it is warm, with good blood circulation. Clean the area with the cleaning cloth. Then allow the cleaning liquid to dry completely.

Step 2

Turn the safety catch of the lancet and remove it. Hold the lancet firmly against the cleaned fingertip and press the purple button. The lancet is triggered and punctures the finger. Gently rub and press the finger towards the punctured area until a sufficiently large drop of blood has formed.

Attention: The lancet can only be triggered once to preserve sterility. Do not test the lancet in advance.

Note: The second lancet included in the kit is a reserve in case the 1st lancet cannot be used. It is not necessary to puncture yourself twice.

Step 3 (Fig. 3, page 33)


Hold the capillary in the center. Make absolutely sure not to close the green piston at the upper end of the capillary. Now hold the capillary from below against the blood drop and wait

until it is automatically drawn into the capillary by capillary force. Fill the capillary up to the white filter with blood (corresponds to 50 μl). To fill the capillary completely with blood, a single drop of blood is not sufficient in most cases. In this case, carefully place the capillary on e.g. the edge of the table to form another blood drop by gently rubbing and pressing your finger and continue filling the capillary.

Note: The capillary will automatically stop drawing up the sample as soon as enough blood has been collected.

Attention: Be careful not to close the small opening at the top, green part of the capillary (e.g. by your finger), otherwise the blood cannot be drawn up by capillary force.

Make sure that the tip of the capillary does not come into contact with anything when you put it down.

No air bubbles must be visible in the capillary, as the correct sample volume is not given. The test cannot be continued in this case.

Step 4

Place the capillary filled up to the white filter e.g. on the edge of the table and apply a plaster to the wound.

Caution: After dressing the wound, proceed to the next step immediately to avoid clotting of the blood sample.


Test procedure

Step 1

Remove the test cassette from the packaging and place it on a flat, dry and light surface.

Attention: Make sure that the test cassette is placed horizontally at all times. **Carry out the next step immediately.**

Sample application

Step 2 (Fig. 2, page 34)

Hold the filled capillary directly over the sample application window (S) of the test cassette. Apply the entire contents of the capillary to the sample application window (S) by slowly pressing down the green piston.


Step 3

Press „Start Timer“ in the **SmarTest Ferritin** app immediately after sample application. A 15-minute timer will start.

Attention: Proceed to the next step immediately.

Step 4

30 seconds after starting this timer, a signal sounds in the app. Now immediately apply 1 hanging drop of chase buffer from a distance of approx. 1 cm to the same sample application window (S) of the test cassette.


Attention: The drop must fall freely from a height of approx. 1 cm to ensure the correct drop volume.

You have 1 minute to apply the chase buffer. Confirm the successful application of the chase buffer in the app. The timer now continues until the 15 minutes of total time have elapsed. A signal sounds after the time has elapsed.

Step 5 (Fig. 5, page 35)

After the timer has expired, you can evaluate the result. Align the template on the display with the test cassette – hold your smartphone parallel to the test cassette. This step must be done within the first 2 minutes after the 15-minute timer has expired in order to obtain a valid result.


Attention: Make sure that the lighting is even. Avoid shadows, light from the side and direct sunlight. **If the time is exceeded, the result is invalid. The measurement must be repeated with a new test set.**

Step 6

The **SmarTest Ferritin** app determines the quantitative result of your measurement based on the ratio of the found intensity of the test band to the control band in comparison with a defined calibration function (4-parameter fit). **This result is displayed directly on your smartphone.** In the **mySmarTest** area, you can display your results at any time and transmit them to healthcare professionals if required. You can configure these settings at any time under „**mySmarTest**“ in the main screen.

Evaluation

The test result is **valid** if:

- a control band (C) appears.
- a control band (C) and a test band (T) appear.

The test result is **invalid** if:

- **no** control band (C) appears.
- **no** control band (C) appears but a test band (T) appears.

Note: If the test result is invalid (missing control band (C)), the test must not be evaluated further. In this case, the test must be performed again with a new test kit. If you have any questions, please contact info@preventis.com.

Interpretation of the test result⁶

Note: The following reference values apply to apparently healthy adults. The reference values for ferritin differ for women and men. Therefore, please orientate towards the biological sex.

Attention: The interpretation of the measured value may differ from this table depending on individual health and physical conditions. The result may be elevated in the presence of infection or inflammatory

events. During pregnancy, there may be changes in threshold values and increased iron requirements. Certain diseases may result in highly elevated ferritin levels as well as altered thresholds. Ferritin levels can rise to ranges of $> 10,000 \mu\text{g/L}$ in extreme cases (e.g. Still's syndrome, severe liver or kidney disease). At such high levels, there is a risk of falsely low test results being displayed. If you have any questions regarding the interpretation of your test result or supplementation, please contact your healthcare professional. If your test result is outside the normal range, seek medical advice for a full evaluation.

	Women	Men
The result is below the normal range. This indicates an iron deficiency	$< 15 \mu\text{g/L}$	$< 15 \mu\text{g/L}$
The result is within the normal range	$15\text{--}200 \mu\text{g/L}$	$15\text{--}300 \mu\text{g/L}$
The result is above the normal range. This indicates an iron oversupply	$> 200 \mu\text{g/L}$	$> 300 \mu\text{g/L}$

Test performance

Analytical and clinical performance

- The test has a measurement range of $10\text{--}500 \mu\text{g/L}$ ferritin in finger prick blood.
- 13 different drugs (e.g. ibuprofen, pantoprazole, aspirin) and food supplements (e.g. 25(OH)-vitamin D_3 , multivitamin tablets, Ferro Sanol) were tested to see if the test reacted to them

(cross-reaction). No relevant cross-reaction could be detected for any of the tested substances. A detailed list is available on request.

- A comparative measurement of 110 patient samples with the reference method turbidimetry (particle-enhanced immunological turbidity test) – a measurement method used in clinical laboratories – showed very good correlation between the two methods (93.6%).
- In studies of accuracy, repeatability and reproducibility, tests were performed with different test batches, smartphones and in different environments. Variations were always smaller than 30%.
- Clinical validity is derived from the agreement between **Preventis SmarTest® Ferritin Home** and an established clinical laboratory test based on turbidimetry, as well as the threshold values outlined in relevant scientific literature.

Limitations of the test

- Although **Preventis SmarTest® Ferritin Home** reliably detects ferritin, false results may occur in individual cases.
- In case of inconclusive results and before implementing therapeutic measures such as supplementation or changing the supplementation dose, it is recommended to contact healthcare professionals.
- As with all diagnostic tests, a conclusive finding should not be made based on a single result, but only after complete clarification by your treating healthcare professional.
- Ferritin levels can rise to ranges of $> 10,000 \mu\text{g/L}$ in extreme cases (e.g. Still's syndrome, severe liver or kidney disease). At such high levels, there is a risk that a falsely low test result will be displayed.

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Literatur/References/Références/Riferimenti/Referanslar

1. **Gao, G. et al. (2019)** Cellular Iron Metabolism and Regulation *Adv Exp Med Biol* 1173: 21–32.
2. **Chiancone, E. et al. (1980)** Ferritin: structural and functional aspects in Radioimmunoassay of Hormones, Proteins and Enzymes, *Proc. Int. Symposium. Excerpta Medica Amsterdam* 197–203.
3. **Knovich, M. A. et al. (2009)** Ferritin for the Clinician *Blood Rev.* 23 (3): 95–104.
4. Iron Factsheet for Health Professionals, National Institutes of Health, Accessed on 11/2021.
5. **Georgieff, M. K. et al. (2020)** Iron deficiency in pregnancy *Am J Obstet Gynecol* 223 (4): 516–524.
6. **Worwood, M. et al. (2016)**: Iron Deficiency Anaemia and Iron Overload, In: *Dacie and Lewis: Practical Haematology*, 12th Ed.

Bei Fragen kontaktieren Sie unseren Kundenservice unter:
If you have any questions, contact our customer service at:
Si vous avez des questions, contactez notre service clientèle à l'adresse suivante :
Se avete domande, contattate il nostro servizio clienti:
Herhangi bir sorunuz varsa, aşağıdaki adresten müşteri hizmetlerimizle iletişime geçin:
+49 6251 70711-0 | info@preventis.com

Immundiagnostik AG
Stubenwald-Allee 8a
64625 Bensheim, Germany

Vertrieb/Distribution/Distribuzione/
Dağıtım:

Preventis GmbH
Stubenwald-Allee 8a, 64625 Bensheim, Germany
T: +49 6251 70711-0 | info@preventis.com
www.preventis.com | shop.preventis.com

2023-09-13