REF PST80100TP (UDI-DI: 04050598008625) | PST80100G5 (UDI-DI: 04050598008632) | PST80100GP (UDI-DI: 04050598008649)

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Intended purpose

Preventis SmarTest Pro® Vitamin D is a rapid test for the quantitative immunological determination of 25(OH)-vitamin D in capillary blood. This value allows monitoring of the vitamin D status and detection of a possible insufficiency or excess. The manual test system consists of an antibody-based strip test, a so-called lateral flow test, including a smartphone app developed for test evaluation. The test is an in vitro diagnostic device and suitable for testing by professional users.

General information

- **Preventis SmarTest Pro® Vitamin D** is intended for in vitro diagnostic use only, which means for use outside of the body.
- Lancets, test cassettes, capillaries and buffer tubes are for single-use only.
- Test components from different production units must not be mixed.
- All test components may only be used for the designated purpose and in combination with the components included in the kit.
- The prepared sample may only be used in combination with the associated test components.
- Sample preparation and test procedure take approximately 5 minutes with additional 15 minutes incubation time. The **SmarTest Vitamin D app** allows parallel incubation of several tests using a multi timer.
- The test must be performed at room temperature.
- Read test instructions carefully before use. Set up the app, take the sample, perform the test and evaluate it exactly according to the test instructions.
- Make sure your smartphone is sufficiently charged (battery level at least 20%).
- The **SmarTest Vitamin D app** requires an active internet connection at regular intervals, but at least every 48 hours, to load the latest calibration data.
- Do not turn off the phone and keep the ringtone on during the application, as the timer uses a signal tone.
- Perform the device qualification and camera test in the same location that the rapid test will be performed in order to accurately simulate lighting conditions. Be sure that the lighting is even. Avoid shadows, light from the sides and direct sunlight.
- Ensure that the camera test card and test cassette are placed on a flat, dry and bright surface, resting horizontally at all times.
- After the 15-minute timer has expired, you have 2 minutes to scan the test cassette with your smartphone. You will only receive a valid result if the test is evaluated within this time period.
- For the security of your data, it will only be transmitted in encrypted form. For exact details please refer to the corresponding Data Privacy Policy from **Preventis GmbH**, which can be found at www.preventis.com.
- Please use only the latest version of your original operating system provided by the manufacturer of your smartphone. A list of supported operating systems can be found at www.link.preventis.com/vitamin-d.

Safety notes

- Do not drink the contents of the buffer tube, avoid direct contact with the liquid. Immediately wash contaminated areas thoroughly after skin contact.
 - If you have problems setting up the app, collecting the specimen, performing or evaluating the test, contact info@preventis.com.
 - Treat the specimen as potentially infectious: do not eat, drink, or smoke during specimen handling. It is recommended to always wash/disinfect hands after testing.
- In case of serious incidents, please contact Preventis GmbH (info@preventis.com) and the appropriate authorities immediately.**

Storage, stability and disposal

- Store the test at room temperature or in the refrigerator (+4 to +30°C), do not freeze.
- The test is sensitive to humidity, direct sunlight and high temperatures. Therefore, protect the test from heat and use immediately after opening the test cassette packaging (within 5 minutes).
- Do not use the test if the packaging or the supplied materials are damaged.
- Do not use the test after the expiration date.
- After use of one test, dispose of all test components except camera test card and test instructions in the residual waste. After using up all the components in the package, please also dispose of the camera test card and test instructions. These must not be used in combination with components from a new package.

Medical background^{1,2,3,4}

„Vitamin D“ refers to a group of molecules that fulfill many different important functions in the body. In contrast to other vitamins, only a small proportion is ingested through food. However, vitamin D can be produced independently by the body when exposed to sunlight. The storage form 25(OH)-vitamin D is measured as a marker for overall vitamin D supply as the majority of total vitamin D consists of this form. Vitamin D is able to bind to almost every cell in the body. Accordingly, vitamin D has various positive effects on the body. It is not only essential for bone formation, but also extremely important for a healthy immune and nervous system, the musculature and the cardiovascular system. Furthermore, vitamin D deficiency has been linked to the development of cancer, diabetes, hypertension, heart attacks and autoimmune diseases. Regular monitoring of vitamin D status is particularly recommended for at-risk groups such as the elderly, people with dark skin color, and those who spend little or no time outdoors.

Test principle

Preventis SmarTest Pro® Vitamin D is an antibody-based strip test for the detection of 25(OH)-vitamin D in capillary blood. The sample is added to a buffer solution and dissolved using a specific sampling system. Then 3–4 drops of the solution are added to the sample application window (S) of the rapid test. In this process, 25(OH)-vitamin D present reacts with gold-labeled anti-25(OH)-vitamin D antibodies and becomes visible as a reddish test band after 15 minutes latest – if the test is performed properly. The appearance of the control band (C) also indicates whether the test has functioned correctly from a technical point of view (see **Evaluation**). Using the **SmarTest Vitamin D app**, the quantitative result is then evaluated based on the color intensity of the control and test bands. The collection and processing of the sample before measurement, as well as the test procedure and evaluation are described step by step in the chapters **Sample extraction**, **Test procedure** and **Evaluation**.

Included materials

PST80100TP (1 test) | PST80100G5 (5 tests) | PST80100GP (20 tests)

- 1 | 5 | 20 **test cassettes** [TEST], individually wrapped with desiccant bag*
- 1 | 5 | 20 **capillaries** [PIP]
- 1 | 5 | 20 **buffer tubes** [BUF]
- 1 | 1 | 1 **camera test card** [CARD]
- 1 | 1 | 1 **instructions for use**
- 1 | 7 | 22 **lancets** [LANCET]

*To prevent the test strip from absorbing moisture, a desiccant bag is included in the packaging. Dispose of the desiccant bag after opening.

Additional materials needed: A smartphone with the **SmarTest Vitamin D app**, an internet connection, plasters, disinfectant.

Setting up the SmarTest Vitamin D app

The **SmarTest Vitamin D app** guides you step-by-step through the sampling and testing process and displays your test result at the end.

Download and installation

Attention: The successful installation (app icon appears on your smartphone) of the **SmarTest Vitamin D app** is mandatory for the evaluation of the rapid test! The **SmarTest Vitamin D app** is available for **iOS**, available in the App Store and **Android**, available in the Google Play Store.

You can find the already supported smartphone models and operating systems at www.link.preventis.com/vitamin-d. If your smartphone is not yet listed, you can qualify it for the test using the camera test card (CARD) (see **Device qualification**).

Please note: Unfortunately, not all smartphones meet the requirements of the app. If you have any questions, please feel free to contact us at info@preventis.com. If you do not receive the device qualification prompt, then go to the next step **Camera test**. Your smartphone model has then already been approved by us.

Device qualification

Attention: Perform the device qualification in the same location that the rapid test will be performed later in order to accurately simulate the lighting conditions. Be sure that the lighting is even. Avoid shadows, light from the sides and direct sunlight.

1. Turn the camera test card (CARD) in such a way that the yellow stripe is on the right side.
 - Note:** If the app does not prompt you to perform the device qualification, go directly to the **Camera test** chapter.
2. Align the template on your display with the six squares on the camera test card. If the alignment is correct, the template will turn green. Make sure to hold your smartphone parallel to the camera test card until you receive a success notification (see Fig. 1a).
 - Note:** Please make sure to hold your smartphone horizontally to the camera test card lying in front of you. Vary the distance of your smartphone to the camera test card until the templates on your screen match the outlines of the squares.


If the device qualification fails repeatedly, check the camera lens for scratches and dirt and remove them if possible. Also check the camera test card for any damage.

Attention: If the device qualification fails repeatedly, your smartphone may not be suitable for evaluation of our rapid tests. For further questions please contact us at info@preventis.com.

Camera test

Attention: Perform the camera test in the same location that the rapid test will be performed later in order to accurately simulate the lighting conditions. Be sure that the lighting is even. Avoid shadows, light from the sides and direct sunlight.

1. Place the camera test card in front of you in such a way that the yellow stripe is facing up and you can read the numbers. Align the template on your display with the image of test cassette ① until you receive a success message (see Fig. 1b).
2. Repeat this step with the image of test cassette ②.

Note: Make sure to hold your smartphone parallel to the camera test card. If the camera test fails repeatedly, check the camera lens for scratches and dirt and remove them if possible. Also check the camera test card for any damage. Please make sure to hold your smartphone horizontally to the camera test card lying in front of you. Vary the distance of your smartphone to the camera test card until the template on your screen matches the outline of the depicted test cassette. If the test continues to fail, your smartphone may not be suitable for evaluation of our rapid tests. For further questions please contact us at info@preventis.com.

Sample extraction


Note: Please allow for approximately 20 minutes to complete the test. Make sure that the test materials are at room temperature. Get the sealed test cassette and all other test components in place. Make sure that the tip of the capillary does not come into contact with anything.

Attention: Leave the test cassette in its packaging.

Check the internet connection of your smartphone and turn on the sound. Open the **SmarTest Vitamin D app** and select „Prepare test“. Then start the animated video instructions for sample extraction, test procedure and evaluation via „To test“.

Please perform all steps one after another without pauses. Do not store the sample in the capillary or in the buffer, but always immediately proceed with the next steps.


1. The lateral areas of the fingertip of the middle or ring finger are particularly suitable for collecting capillary blood. Rub the selected finger until it is warm, with good blood circulation. Disinfect the area, then allow the disinfectant to dry completely.
2. Turn the safety catch of the lancet and remove it. Hold the lancet firmly against the cleaned fingertip and press the yellow button. The lancet is triggered and punctures the finger. Gently rub and press the finger towards the punctured area until a sufficiently large drop of blood has formed.
 - Attention:** The lancet can only be triggered once to preserve sterility. Do not test the lancet in advance.
3. Hold the capillary horizontally or from below against the blood drop and wait until the blood drop is automatically drawn into the capillary by capillary force (see Fig. 2a). Fill the capillary completely with blood (corresponds to 10 µl) without removing it from the finger.
 - Attention:** No air bubbles must be visible in the capillary, as the correct sample volume is not given. The test cannot be continued if this happens.
 - Please proceed to the next step immediately** to avoid clotting of the blood sample.
4. Place the buffer tube on a flat, dry surface with the foil facing upwards. Place the filled capillary to the foil of the buffer tube at a slight angle and pierce the foil carefully. Press the capillary and the buffer tube firmly together until the notch of the capillary is completely resting on the edge of the buffer tube (see Fig. 2b).
 - Note:** If the capillary comes off while being pressed together with the buffer tube, please put it back on.
 - Attention:** Please proceed immediately to the next step to avoid clotting of the blood sample.
5. Shake the buffer tube with the capillary attached until blood and buffer are well mixed.
6. Now you can apply a plaster to the punctured area.


Test procedure

1. Remove the test cassette from the packaging and place it on a flat, dry and light surface (Fig. 3).

Attention: Make sure that the test cassette is placed horizontally at all times. **Carry out the next step immediately.**

2. Remove the round cap from the buffer tube/capillary with a slight twisting motion and hold the tube approx. 1 cm above the sample application window (S). Drop 3–4 hanging drops onto the sample application window (S) of the test cassette by squeezing the center of the tube (siehe Fig. 4a).
 - Attention:** The drops must fall freely from a height of approx. 1 cm to ensure the correct drop volume.
3. Press „Start Timer“ in the **SmarTest Vitamin D app** immediately after sample application. A 15-minute timer starts. After expiration, a signal sounds.
 - Attention:** If you want to perform several tests at once, please choose the multi-timer function on the welcome screen.
4. After the timer has expired, you can evaluate the result. Align the template on the display with the test cassette – hold your smartphone parallel to the test cassette. This step must be done within the first 2 minutes after the 15-minute timer has expired in order to obtain a valid result (see Fig. 4b).
 - Attention:** Make sure that the lighting is even. Avoid shadows, light from the side and direct sunlight. **If the time is exceeded, the result is invalid. The measurement must be repeated with a new test set.**
5. The **SmarTest Vitamin D app** determines the quantitative result of your measurement based on the ratio of the found intensity of the test band to the control band in comparison with a defined calibration function (4-parameter fit). **This result is displayed directly on your smartphone.**

Evaluation

The test result is **valid** if:

- a control band (C) appears
- a control band (C) and a test band (T) appear.

The test result is **invalid** if:

- no control band (C) appears.
- no control band (C) appears but a test band (T) appears.

Note: If the test result is invalid (missing control band (C)), the test must not be evaluated further. In this case, the test must be performed again with a new test kit. If you have any questions, please contact info@preventis.com.

Interpretation of the test result⁵

< 30 ng/ml	The result indicates insufficient supply
30–100 ng/ml	The result indicates an optimal supply
> 100 ng/ml	The result indicates an oversupply

Attention: The interpretation of the test result may differ from this table depending on individual health and physical conditions. In individual cases, the intake of high-dose supplements may lead to falsely low test results.

Test performance

Analytical performance

- The test detects 25(OH)-vitamin D₂, 25(OH)-vitamin D₃, 1,25(OH)-vitamin D₂, 1,25(OH)-vitamin D₃, 3-epi-25-OH-vitamin D₂, 3-epi-25-OH-vitamin D₃, and 24R,25-OH-vitamin D₃. Due to the normal (relatively low) abundance of these forms of vitamin D the impact is clinically not relevant.
- The test has a measurement range of 5–100 ng/mL.
- The **Preventis SmarTest Pro® Vitamin D** test result is not affected by typical substances found in blood (e.g. cholesterol), dietary supplements (e.g. vitamin B₁₂ or vitamin C), or blood-thinning medications (Aspirin, Xarelto). A detailed list of all tested substances is available upon request from **Preventis GmbH**.
- A comparative measurement of 27 patient samples with the reference method LC-MS/MS (liquid chromatography-mass spectrometry/mass spectrometry) – a measurement method used in clinical laboratories – showed a very good correlation between the two methods (87.5%).
- For repeated measurements with one rapid test batch, the coefficient of variation is 25.8%, for repeated measurements with different rapid test batches it is 3–37%.

Limitations of the test

- Although **Preventis SmarTest Pro® Vitamin D** reliably detects 25(OH)-vitamin D, false results may occur in individual cases.
- As with all diagnostic tests, a conclusive finding should not be made based on a single result.
- In individual cases, the intake of high-dose supplements may lead to falsely low test results.

References

1. Nair R, Maseeh A (2012) Vitamin D: The “sunshine” vitamin. J Pharmacol Pharmacother. 3(2): 118–126
2. Holick MF (2012) Vitamin D: extraskeletal health. Rheum Dis Clin North Am. 38(1): 141–60
3. Wen H, Baker JF (2011) Vitamin D, immunoregulation, and rheumatoid arthritis. J Clin Rheumatol. 17(2): 102–7
4. Lips P, et al. (2019) Current vitamin D status in European and Middle East countries and strategies to prevent vitamin D deficiency: a position statement of the European Calcified Tissue Society Eur J Endocrinol. 180(4): P23–P54
5. Holick et al. (2011) Evaluation, Treatment, and Prevention of Vitamin D Deficiency: an Endocrine Society Clinical Practice Guideline J Clin Endocrinol Metab. 96(7): 1911–30

Short Instructions Preventis SmarTest Pro® Vitamin D

1. Disinfect the fingertip, then allow the disinfectant to dry completely.
2. Puncture the fingertip using the lancet.
3. Fill the capillary completely with blood.
4. Insert capillary into buffer tube, press firmly together and shake.
5. Remove the test cassette from the packaging.
6. Remove the cap from buffer tube/capillary. Drop 3–4 hanging drops onto the sample application window (S) of the test cassette.
7. Allow 15 minutes incubation time.
8. After the signal tone, evaluate the result using the **SmarTest Vitamin D app** within 2 minutes.
9. The result is displayed directly on the smartphone.

CE [ND] Rev. 2022-05-24

- Temperature limit
- Manufacturer
- To be used with
- Batch number
- Catalogue number
- Please see reverse
- Use by
- Follow test instructions
- Contents sufficient for <n> tests
- Not reusable
- In vitro diagnostic device: only to be used outside of the body
- Sterilized using irradiation
- Keep away from sunlight

PREVENTIS

Distribution:

Preventis GmbH
 Stubenwald-Allee 8a
 64625 Bensheim, Germany
 T: +49 6251 70711-0
 F: +49 6251 70711-299
info@preventis.com
www.preventis.com
www.shop.preventis.com

Immundiagnostik AG
 Stubenwald-Allee 8a
 64625 Bensheim, Germany