→REF KST15110TP
KST15110G3
KST15110GP

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→REF KST15110G3
KST15110GP

Intended Use

PreventID® Bladder is an immunochromatographic rapid test for the parallel determination of the bladder cancer marker (BCM, human complement factor H-related protein) and hemoglobin in urine for an early detection of bladder cancer.

Test Principle

The PreventID® Bladder test employs a unique combination of monoclonal and polyclonal antibodies to selectively identify BCM – a specific bladder cancer protein – in urine samples. BCM is a well recognized prognostic marker for the detection of primary and recurrent bladder tumors.

Intermittent hematuria is a common symptom of bladder cancer. The PreventID® Bladder employs specific antibodies to determine even smallest amounts of hemoglobin (≥ 250 ng/mL) in urine samples.

The parallel detection of BCM and hemoglobin improves the rate of cancer detection significantly.

Materials Provided

Materials Provided:

- test devices (with sample droppers), individually packed **TEST**
- manual

Materials Required but not Provided: Timer or stop watch

Storage and Stability

The PreventID® Bladder test should be stored at room temperature (2–30°C). The test device is sensitive to humidity and heat. Therefore, perform the test immediately after removing the test device from the pouch. Do not use beyond the expiry date.

Precautions

1. For *in vitro* diagnostic use only.
2. Do not eat or smoke while handling specimen.
3. Wear protective gloves and wash hands thoroughly after performing the test.
4. Avoid splashing or aerosol formation while handling specimen and performing the test.
5. All samples and materials used should be treated as potentially infectious and disposed in a biohazard container. Clean all contaminated objects and surfaces carefully.
6. Do not use test if the pouch is torn or if the membrane of the test device is visibly damaged. Note the expiry date.
7. Read the instruction carefully before performing the test.
8. Do not mix reagents from different lots.
9. If you have any questions please contact Preventis GmbH.

Specimen Collection

Specimen collection should not be performed during or within three days of a menstrual period. Subjects who are on aspirin medication and obtain a positive result should repeat the test after aspirin has been stopped for 7 days.

1. Only urine specimens should be used in this assay. Avoid sample contamination with toilet water.
2. Specimens that are not tested the same day can be refrigerated and stored at 4 to 8°C up to 3 days.

Test Procedure

1. Remove the test device from the pouch and place it on a flat, dry surface.
2. If the specimen has been refrigerated, bring it to room temperature.
3. Transfer 2 to 3 drops of urine to each sample application window of the test device (Fig. 1).
4. As the test panel begins to work, you will see purple colour move across the respective result window in the center of the test device.
5. Interpret test results at 9 minutes. Do not interpret test results after 10 minutes.


Fig. 1: PreventID® Bladder Test device

Caution: The above interpreting time is based on reading the test results at room temperature of 15 to 30°C. If your room temperature is significantly lower than 15°C, then the interpreting time should be properly increased.

References

1. Badrinath R. Konety and Robert H. Getzenberg: Urine based markers of urological malignancy. The Journal of Urology, Vol. 165, 600-611, Feb, 2001
2. MP Raitanen, E Kaasinen, E Rintala, E Hansson, P Nieminen, R Aine, TLJ Tammela and The Finn Bladder Group: Prognostic utility of human complement factor H related protein test. British Journal of Cancer (2001) 85(4), 552-556

Test Interpretation (Fig. 2)

A colour line will appear at the left section of the result window to show that the test is working properly. This line is the control line „C“. A colour test line will appear at „T“ of the Hb section of the result window (indicated by the letters „Hb“) when the hemoglobin level in a sample is 250 ng/mL or higher. A colour test line will appear at „T“ of the Bladder Cancer Marker section of the result window (indicated by the letters „BCM“) when the bladder cancer marker is detected in the urine.

Positive: The presence of the „C“ colour line along with either or both visible test lines „T“ of the Hb and BCM sections (no matter which line appears first) indicates a positive result (Fig. 2a, 2b, 2c).

Negative: The presence of only the „C“ colour line within the result windows in both the Hb and BCM sections indicates a negative result (Fig. 2d).

Invalid: If after performing the test no colour line or no „C“ colour line is visible within the result window the test is considered invalid. It is recommended that the specimen be re-tested (Fig. 2e, 2f).


Fig. 2: PreventID® Bladder Test Interpretation

A positive result in both Hb and BCM indicates a likelihood of bladder cancer, however requires further confirmation.

A positive result only in Hb: Other common causes of blood in urine need to be ruled out, such as bleeding from stones or enlarged prostate, and the subject needs to be followed up with repeated testing.

Negative results do not exclude bladder cancer, other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

Note: Once a positive result has been established (at 9 minutes), the result will not change. However, in order to prevent any incorrect results, the test result should not be interpreted after 10 minutes.

Clinical Evaluation

Clinical study

A total of 225 patients with biopsy confirmed bladder cancer took part in this study, of the 225 patients, stages were determined on 223, and grades are determined on 210 patients.

Stage	N	Sensitivity (BCM only)	Sensitivität (BCM and Hb)
Ta	14	57.0%	57.0%
T1	69	82.6%	85.5%
T2, T3	140	86.4%	94.2%
Overall	223	83.4%	89.2%
Grade	N		
1	47	38.2%	42.5%
2	61	63.9%	70.5%
3	102	89.2%	97.1%
Overall	210	70.5%	77.0%

PreventID® Bladder had overall sensitivity varying from 77% (WHO Grading) to 89.2% (TNM Staging).

Relative Specificity

To obtain Clinical Specificity data, urine from 38 healthy individuals with no known bladder tumors are tested:

N	True Negative
Bladder/hb Positive	4
Bladder/hb Negative	34
Total	38

PreventID® Bladder had overall specificity of 89.5%. In an in house study, no false positive was observed with microalbumin up to 40 µg/mL, however, this is not clinically verified.

Test Limitations

Test results are only reliable if you follow the instructions for use carefully. The test is limited to the detection of bladder cancer in human urine. Use the test only once. Although the test is very accurate, a low incidence of false results can occur. If negative or questionable results are obtained, the test should be repeated on a fresh urine specimen using a new test device.

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

US: all products: Research Use Only. Not for use in diagnostic procedures.

Temperature limitation	Manufacturer	CE	Status: 2020-08-25
In vitro diagnostic device	Lot number		
Catalogue number	Expiry date	Distributed by: Preventis GmbH Stubenwald-Allee 8a 64625 Bensheim, Germany Phone: +49 6251 70711-0 Fax: +49 6251 70711-299 info@preventis.com www.preventis.com	
To be used with	Do not reuse		
Read user instructions	Contains sufficient for <n> tests		